Supplementary Material

Supplementary Figures



Figure S1. Comparison of TrkB-T1 protein levels in primary cultures of mixed cortical cells or astrocytes. (A) Representative immunoblot analysis of equivalent amounts of protein lysates corresponding to cultures enriched, respectively, in neurons (Neu) or astrocytes (Astro). Levels of the TrkB-T1 isoform and β -actin, used for normalization, are shown. (B) Quantitation by densitometric analysis of normalized TrkB-T1 levels showing the fold-increase in astrocytes compared to the mixed cultures. Means \pm SEM are presented and differences were analyzed by Student's t-test (n = 3)



Figure S2. Biological processes associated with $Bio-sTT1_{Ct}$ protein interactions. Representation of the most relevant Gene Ontology (GO) terms related to biological processes for selected proteins interacting with $Bio-sTT1_{Ct}$ in basal conditions (A) or after treatment with BDNF (B) or NMDA (C). The dot size represents the number of proteins from our data set related to each process. Dots are colored according to their statistical significance, which is set by a color scale referring to -log10 (adjusted p-value).



Figure S3. Total RhoGDI1 levels are not affected by BDNF or NMDA treatment. (A) Western blot analysis of cortical cultures incubated with Bio-sTT1_{Ct} or Bio-TMyc (25 μ M) for 30 min before treatment with BDNF (100 ng/ml) or NMDA for 30 min. A representative experiment is presented showing RhoGDI levels in total lysates. (B) Quantitation by densitometric analysis of RhoGDI total levels. Means ± SEM are presented relative to basal conditions with Bio-TMyc (100%). Data were analyzed by ANOVA test followed by *posthoc* Tukey's test, n = 6-7.



Figure S4. Neuroprotection due to $TT1_{Ct}$ preincubation is maintained when applied at the time of damage induction. (A) MTT assay design to evaluate $TT1_{Ct}$ neuroprotection in primary cortical cultures treated with CPPs at different time points of damage induction. TMyc or $TT1_{Ct}$ were added 30' before (CPPs pre), at the same time (CPPs 0), 30 min (CPPs 30) or 60 min (CPPs60) after treatment with NMDA (100 μ M). Neuronal viability was established after 4 h of excitotoxicity. (B) Neuronal viability in cultures treated as above indicated. Individual results and means \pm SEM are presented relative to the values obtained for the untreated cells (100%). Data were analyzed using two-way ANOVA test followed by *post hoc* Bonferroni test, n = 4-8.



Figure S5. Bio-TT1_{Ct} and Bio-TMyc can cross the BBB of undamaged animals. Analysis of biotinylated-peptide delivery to undamaged male mice cortex. Bio-TT1_{Ct} and Bio-TMyc (3 nmol/g) were injected, and animals were sacrificed 30 min after. Detection was made in coronal sections by Fluorescein Avidin D (green). Peptide delivery was observed in cell projections (arrowheads) and bodies (asterisks) of cortical cells. Representative confocal microscopy images of cortical areas correspond to single sections. Scale bar, 20 μ m.



Figure S6. Analysis of the expression of TrkB isoforms and cell-specific markers in damaged brain. Immunohistochemistry of brain coronal sections prepared from male animals i.v. injected with Bio-TMyc (3 nmol/g) (A) or vehicle (B) one hour after damage initiation and sacrificed after 5 h of injury. Sections were stained with isoform-specific TrkB antibodies (TrkB-FL and TrkB-T1), NeuN, GFAP and DAPI. Representative images of cortical areas of the infarct border, obtained with a Cell Observer, are shown. Arrows indicate concurrent increased expression of TrkB-T1 and GFAP in astrocytes. Scale bar, 50 µm.



Figure S7. Leakage of mouse immunoglobulins to the brain cortex due to BBB breakage is decreased in $TT1_{Ct}$ -treated animals. Brain coronal sections of male animals sacrificed 24 h after insult were analyzed by immunohistochemistry with an anti-mouse secondary antibody. Heavy staining of blood vessels, infiltrated cells and high backgrounds were observed in the ischemic brain of animals injected with the control peptide and decreased in $TT1_{Ct}$ -treated animals. Representative images correspond to single sections. Scale bar, 20 µm.

Table S1. Key reagents and resources with indication of sources and identifiers.

Antibodies Mouse monoclonal anti-β-Actin Sigma-Aldrich Cat#A5441 RRID: AB 10920058 Rabbit polyclonal anti-C3d Dako Cat#A0636 Mouse monoclonal anti-C068 (ED1) Millipore Cat#MA535: RRID: AB 177576 Rabbit polyclonal anti-GFAP Millipore Cat#G6519; RRID: AB 10153 Rabbit monoclonal anti-GFAP Millipore Cat#G6171 Raubit monoclonal anti-GFAP Millipore Cat#G6171 Runo AB 11212597 Mouse monoclonal anti-GFAP Millipore Cat#G6171 Runo AB 14452 RRID: AB 39504 Cat#G6171 Runo AB 14452 BD Biosciences Cat#G6171 Runo AB 14452 BD Biosciences Cat#G6171 Runo AB 14452 BD Biosciences Cat#G1741 Runo AB 14455 RRID: AB 39504 RRID: AB 39504 Mouse monoclonal anti-RBFOX3/NeuN Novus Cat#G1741 Rabbit polyclonal anti-RBFOX3/NeuN Novus Cat#G1741 Rabbit polyclonal anti-RBFOX3/NeuN Novus Cat#G17419269358; Rabbit polyclonal anti-RBFOX3/NeuN Notuse Cat#G17421263263 Rabbit poly	REAGENT or RESOURCE	SOURCE	IDENTIFIER
Mouse monoclonal anti-β-Actin Sigma-Aldrich Cat#/A541 RRID: AB 10920058 Rabbit polyclonal anti-C3d Dako Cat#A0063 RRID: AB 578478 Mouse monoclonal anti-CD68 (ED1) Millipore Cat#0A5175 Rabbit polyclonal anti-CD68 (ED1) Millipore Cat#0A515; RRID: AB 17576 Rabbit polyclonal anti-GFAP Millipore Cat#06-519; RRID: AB 11212597 Mouse monoclonal anti-GFAP Millipore Cat#06171 Rubit polyclonal anti-Iba1 Wako Cat#019-19741 Mouse monoclonal anti-MEF2D BD Biosciences Cat#019-19741 Mouse monoclonal anti-RBFOX3/NcuN Novus Cat#019-19741 Mabit polyclonal anti-RBFOX3/NcuN Novus Cat#A0511; RRID: AB 398095 Rabbit polyclonal anti-RBFOX3/NcuN Novus Cat#A051; RRID: AB 1103747 Rabbit polyclonal anti-TrkB-FL (C-ter region) Santa Cruz Cat#A261; RRID: AB 632554 Rabbit polyclonal anti-TrkB-T1 (isoform-specific C-ter) Custor made Cat#A2-11236 Goat anti-mouse IgG Alexa Fluor 546 Molecular Probes Cat#A-11030; RRID: AB 223505 Goat anti-rabbit IgG Alexa Fluor 546 Molecular Probes Cat#A-21236 <td< th=""><th>Antibodies</th><th></th><th></th></td<>	Antibodies		
RRID: AB 10920058 Rabbit polyclonal anti-C3d Dako Cat#A0063 RID: AB 578478 Mouse monoclonal anti-C568 (ED1) Millipore Cat#MA1435; RRID: AB 177576 Rabbit polyclonal anti-pS133-CREB Millipore Cat#MA153; RRID: AB 1121257 Rabbit monoclonal anti-GFAP Millipore Cat#G6171 RRID: AB 1121257 RRID: AB 1840893 Rabbit polyclonal anti-GFAP Millipore Cat#G6171 RRID: AB 1840893 RRID: AB 1840893 Rabbit polyclonal anti-GFAP Millipore Cat#G171 RRID: AB 1840893 RRID: AB 1840893 RRID: AB 1840893 Mouse monoclonal anti-MEF2D BD Biosciences RRID: AB 1840893 Mouse monoclonal anti-RBFOX3/NeuN Novus Cat#AD951; RRID: AB 22930 Rabbit polyclonal anti-reuronal-specific enolase (NSE) Millipore Cat#AB 22751 Rabbit polyclonal anti-TrkB-FL (C-ter region) Santa Cruz Cat#AA 20120; RRID: AB 22751 Rabbit polyclonal anti-rabbit 12G Alexa Fluor 546 Molecular Probes Cat#A-11030; RRID: AB 235805 Goat anti-rabbit 12G Alexa Fluor 546 Molecular Probes Cat#A1-1030; RRID: AB 235805 Goat	Mouse monoclonal anti- β -Actin	Sigma-Aldrich	Cat#A5441
Rabbit polyclonal anti-C3d Dako Cat#MAB1435; RRID: AB 578478 Mouse monoclonal anti-CD68 (ED1) Millipore Cat#MAB1435; RRID: AB 578478 Rabbit polyclonal anti-pS133-CREB Millipore Cat#MAB1435; RRID: AB 177576 Rabbit polyclonal anti-GFAP Millipore Cat#MAB3600 RRID: AB 11212597 RRID: AB 11212597 RRID: AB 11212597 Mouse monoclonal anti-GFAP Millipore Cat#MAB3600 RRID: AB 839504 RRID: AB 839504 Mouse monoclonal anti-Iba1 Wako Cat#O19-19741 RRID: AB 839504 RRID: AB 39504 Mouse monoclonal anti-REF2D BD Biosciences Cat#A18971-9269358; RuD: AB 198095 RRID: AB 1980951; RRID: AB 1980951; Rabbit polyclonal anti-neuronal-specific enolase (NSE) Millipore Cat#A18971; Rabbit polyclonal anti-REF2D Santa Cruz Cat#A-2136 Rabbit polyclonal anti-Rer region of RhoGDI Santa Cruz Cat#A-11030; Rabbit polyclonal anti-TrkB-FL (C-ter region) Santa Cruz Cat#A-11030; RuD: AB 144695 Goat anti-mouse IgG Alexa Fluor 546 Molecular Probes Cat#A-11030; Goat anti-rabbit IgG Alexa Fluor 546 Molecular Probes Cat#A11035 Goat anti-rabbit IgG Alexa Fluor 546 Molecular Probes Cat#A11034 RRID: AB 14	,	0	RRID: AB 10920058
Relit	Rabbit polyclonal anti-C3d	Dako	Cat#A0063
Mouse monoclonal anti-CD68 (ED1) Millipore Catt/MAB1435; RABbit polyclonal anti-pS133-CREB Millipore Catt/MAB1435; RABbit monoclonal anti-GFAP Rabbit monoclonal anti-GFAP Millipore Catt/MAB360 RRID:AB 11212597 Mouse monoclonal anti-GFAP Millipore Catt/MAB360 RRID:AB 11212597 Mouse monoclonal anti-GFAP Millipore Catt/MAB360 RRID:AB 1840893 Rabbit polyclonal anti-Iba1 Wako Catt/MOB19-19741 RRID:AB 1840893 Mouse monoclonal anti-MEF2D BD Biosciences Catt/MAB 398095 Mouse monoclonal anti-RBFOX3/NeuN Novus Catt/MBP1-9269385; Mabbit polyclonal anti-neuronal-specific enolase (NSE) Millipore Catt/AB951; Rabbit polyclonal anti-TrkB-FL (C-ter region) Santa Cruz Catt#AB91; Rabbit polyclonal anti-TrkB-FL (C-ter region) Santa Cruz RRID:AB 632554 Rabbit polyclonal anti-TrkB-FL (G-ter region) Santa Cruz Catt#A-11030; RRID:AB 103747 Molecular Probes Catt#A-11030; RRID:AB 227516 Molecular Probes Catt#A-11030; Goat anti-mouse IgG Alexa Fluor 546 Molecular Probes Catt#A-11030;			RRID: AB_578478
Rabbit polyclonal anti-pS133-CREB Millipore RtID: AB 177576 Rabbit monoclonal anti-GFAP Millipore Cat#06-519; RRID: AB 310153 Rabbit monoclonal anti-GFAP Millipore Cat#MAB360 Mouse monoclonal anti-Iba1 Wako Cat#019-19741 Mouse monoclonal anti-MEF2D BD Biosciences Cat# 610774, RRID: AB 398095 Mouse monoclonal anti-MEF2D BD Biosciences Cat# 610774, RRID: AB 398095 Mouse monoclonal anti-RBFOX3/NeuN Novus Cat#NPB1-9269385; RRID: AB 398095 Mouse monoclonal anti-RBFOX3/NeuN Novus Cat#NPB1-9269385; RRID: AB 92390 Rabbit polyclonal anti-neuronal-specific enolase (NSE) Millipore Cat#NPB1-9269385; RRID: AB 922701 Rabbit polyclonal anti-TrkB-FL (C-ter region) Santa Cruz Cat#Se-11; RRID: AB 2227516 Rabbit polyclonal anti-TrkB-FL (G-ter region) Santa Cruz Cat#A-11030; RRID: AB 2355805 Goat anti-mouse IgG Alexa Fluor 546 Molecular Probes Cat#A-11030; RRID: AB 2535805 Goat anti-rabbit IgG Alexa Fluor 546 Molecular Probes Cat#A-11035 Goat anti-rabbit IgG-heavy and light chain-HRP Bethyl Cat#A1035 RRID: AB 10892625 Donkey anti-	Mouse monoclonal anti-CD68 (ED1)	Millipore	Cat#MAB1435;
Rabbit polyclonal anti-pS133-CREB Millipore Catt/06-519; RID:AB 310153 Rabbit monoclonal anti-GFAP Millipore Catt/06-519; RID:AB 11212597 Mouse monoclonal anti-GFAP Millipore Catt/06-519; RID:AB 11212597 Mouse monoclonal anti-Iba1 Wako Catt/06-6171 RRD:AB 1840893 Mouse monoclonal anti-Iba1 Wako Catt/0174, RRD:AB 398095 Mouse monoclonal anti-REF2D BD Biosciences Catt/0174, RRD:AB 398095 Mouse monoclonal anti-RBFOX3/NeuN Novus Catt/0774, RRD:AB 3980951; RRD:AB 103747 Rabbit polyclonal anti-neuronal-specific enolase (NSE) Millipore Catt/AB951; RRD:AB 92390 Rabbit polyclonal anti-TerRefice C-ter region) Santa Cruz Catt#sc-11; RRD:AB 92390 Rabbit polyclonal anti-TrkB-FL (C-ter region) Santa Cruz Catt#sc-11; RRD:AB 22554 Rabbit polyclonal anti-TrkB-FL (C-ter region) Santa Cruz Catt#sc-11; RRD:AB 23554 Goat anti-mouse IgG Alexa Fluor 546 Molecular Probes Catt#A-11030; RRD:AB 2535805 Goat anti-mouse IgG Alexa Fluor 647 Molecular Probes Catt# A11034 RRD:AB 2535805 Goat anti-rabbit IgG Alexa Fluor 546 Molecular Probes Catt# A10034 RRD:AB 2535805 Goat anti-rabbit IgG Alexa Fluor 546 Molecular Probe			RRID: AB 177576
Rabbit monoclonal ani-GFAP Millipore Cat#MAB360 RRID:AB 11212597 Mouse monoclonal anti-GFAP Millipore Cat#A6171 RRID:AB 11212597 Mouse monoclonal anti-GFAP Millipore Cat#A6171 RRID:AB 840893 Mouse monoclonal anti-Ibal Wako Cat#A6171 RRID:AB 398095 Mouse monoclonal anti-MEF2D BD Biosciences Cat#A619-19741 RRID:AB 398095 Mouse monoclonal anti-RBFOX3/NcuN Novus Cat#A81951; RRID:AB 398095 Mouse monoclonal anti-reuronal-specific enolase (NSE) Millipore Cat#A8951; RRID:AB 92390 Rabbit polyclonal anti-TrkB-FL (C-ter region) Santa Cruz Cat#A8951; RRID:AB 92207516 Rabbit polyclonal anti-TrkB-TL (isoform-specific C-ter) Custom made Cat#A-11030; RRID:AB 143656 Goat anti-mouse IgG Alexa Fluor 546 Molecular Probes Cat#A-11030; RRID:AB 2535805 Goat anti-rabbit IgG Alexa Fluor 546 Molecular Probes Cat#A-11034 RRID:AB 2535217 Goat anti-rabbit IgG Alexa Fluor 546 Molecular Probes Cat#A11035 RRID:AB 2536217 Goat anti-rabbit IgG Alexa Fluor 546 Molecular Probes Cat#A11035 RRID:AB 143051 Donkey anti-rabbit IgG Alexa Fluor 546 Molecular Probes Cat#A11035 RRID:AB 143051 Donkey anti-rabbit IgG-heavy and light chain-HRP	Rabbit polyclonal anti- pS133-CREB	Millipore	Cat#06-519;
Rabbit monecional ani-GrAPMilliporeCat#MAB360Mouse monecional anti-GFAPMilliporeCat#G6171Rabbit polyclonal anti-Iba1WakoCat#G1671Mouse monecional anti-MEF2DBD BiosciencesCat#G1674,Mouse monecional anti-MEF2DBD BiosciencesCat#G1674,Rabbit polyclonal anti-neuronal-specific enolase (NSE)MilliporeCat#G1674,Rabbit polyclonal anti-reuronal-specific enolase (NSE)MilliporeCat#A8951;Rabbit polyclonal anti-TrkB-FL (C-ter region)Santa CruzCat#A8951;Rabbit polyclonal anti-TrkB-FL (C-ter region)Santa CruzCat#A6254Rabbit polyclonal anti-TrkB-FL (C-ter region)Santa CruzCat#A64554Goat anti-mouse IgG Alexa Fluor 546Molecular ProbesCat#A-11030;Goat anti-mouse IgG Alexa Fluor 546Molecular ProbesCat#A-11030;Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat#A11035Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat#A11035Donkey anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat#A11035Chemical and biological productsAra C (usda 110 µM)Sigma-AldrichCat#A1008;RHD:AB 143051Conkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A120-108P;RHD:AB 121256Cat#A1103RRD:AB 143051Cat#A100-137P;RHD:AB 100 ng/ml)Sigma-AldrichCat#C1768;B27 serum free supplementThermo FisherCat#17504044Recombinant human/murinc/rat brain-derived neurotrophic Icat region floo ong/ml)Cat#A100311		N ('11'	<u>RRID:AB_310153</u>
Mouse monoclonal anti-GFAP Millipore Cath/G6171 Rabbit polyclonal anti-Iba1 Wako Cat#019-19741 Rabbit polyclonal anti-Iba1 Wako Cat#019-19741 Mouse monoclonal anti-MEF2D BD Biosciences Cat# 610774, RRID-AB 398095 RRID-AB 398095 RRID-AB 398095 Mouse monoclonal anti-RBFOX3/NeuN Novus Cat#AB951; RRID:AB 22390 RRID:AB 22390 RRID:AB 22390 Rabbit polyclonal anti-neuronal-specific enolase (NSE) Millipore Cat#AB951; RRID:AB 22320 RRID:AB 2227516 RRID:AB 2227516 Rabbit polyclonal anti-TrkB-FL (C-ter region) Santa Cruz Cat#A-1103 (RRID:AB 632554 Rabbit polyclonal anti-TrkB-FL (C-ter region) Santa Cruz Cat#A-11030; Ratin-mouse IgG Alexa Fluor 546 Molecular Probes Cat#A-11030; Goat anti-mouse IgG Alexa Fluor 546 Molecular Probes Cat#A11034 Goat anti-rabbit IgG Alexa Fluor 546 Molecular Probes Cat# A11034 RRID: AB 12376217 Goat anti-rabbit IgG Alexa Fluor 546 Molecular Probes Cat#A11034 Donkey anti-rabbit IgG-heavy and light chain-HRP Bet	Rabbit monocional ani-GFAP	Millipore	Cat#MAB360
Monse monoclonal anti-GrAt Minipole RRID:AB 1840893 Rabbit polyclonal anti-Iba1 Wako Cat#019-19741 Mouse monoclonal anti-MEF2D BD Biosciences Cat# 610774, RRID:AB 398095 Mouse monoclonal anti-RBFOX3/NeuN Novus Cat#AB1103747 Rabbit polyclonal anti-neuronal-specific enolase (NSE) Millipore Cat#AB951; RRID:AB 20390 Rabbit polyclonal anti-TrkB-FL (C-ter region) Santa Cruz Cat#Acs-360 RRID:AB 632554 Rabbit polyclonal anti-TrkB-FL (C-ter region) Santa Cruz Cat#Ac-11030; RRID:AB 632554 Rabbit polyclonal anti-TrkB-FL (isoform-specific C-ter) Custom made Cat#A-11030; RRID:AB 2335805 Goat anti-mouse IgG Alexa Fluor 546 Molecular Probes Cat#A-11034 RRID:AB 2353805 Goat anti-rabbit IgG Alexa Fluor 647 Molecular Probes Cat# A11034 RRID:AB 2535805 Goat anti-rabbit IgG Alexa Fluor 546 Molecular Probes Cat# A11034 RRID:AB 12876217 Goat anti-rabbit IgG-heavy and light chain-HRP Bethyl Cat#A11034 RRID:AB 12892625 Donkey anti-rabbit IgG-heavy and light chain-HRP Bethyl Cat#A11034 RRID:AB 10892625 Donkey anti-rabbit IgG-heavy and light chain-HRP Bethyl Cat#A11044 Recombinant human/murie/rat brain-derived neurotrophic fact	Mouse monoclonal anti CEAP	Millipore	Cat#G6171
Rabbit polyclonal anti-Iba1 Wako Cat#019-19741 Mouse monoclonal anti-MEF2D BD Biosciences Cat# 610774, RRID: AB 393004 Mouse monoclonal anti-MEF2D BD Biosciences Cat# 610774, RRID: AB 398095 Mouse monoclonal anti-MEF2D Novus Cat#NBP1-92693SS; RRID: AB 1103747 Rabbit polyclonal anti-neuronal-specific enolase (NSE) Millipore Cat#AB951; RRID: AB 2227516 Rabbit polyclonal anti-TrkB-FL (C-ter region) Santa Cruz Cat#sc-11; RRID: AB 2227516 Rabbit polyclonal anti-TrkB-FL (C-ter region) Santa Cruz Cat#sc-11; RRID: AB 632554 Rabbit polyclonal anti-TrkB-FL (C-ter region) Santa Cruz Cat#sc-11; RRID: AB 632554 Goat anti-mouse IgG Alexa Fluor 546 Molecular Probes Cat#A-11030; RRID: AB 2535805 Goat anti-rabbit IgG Alexa Fluor 647 Molecular Probes Cat# A11034 RRID: AB 2576217 Goat anti-rabbit IgG Alexa Fluor 546 Molecular Probes Cat# A1034 RRID: AB 2576217 Goat anti-rabbit IgG-heavy and light chain-HRP Bethyl Cat# A11035 RRID: AB 1292625 Donkey anti-rabbit IgG-heavy and light chain-HRP Bethyl Cat# A11037 RRID: AB 1292625 Donkey anti-rabbit IgG-heavy and light chain-HRP Bethyl Cat# A11034 RRID: AB 1292625 Donkey anti	Mouse monocional anti-OFAI	winipore	RRID: AR 1840893
IndexRRID: AB 39504Mouse monoclonal anti-MEF2DBD BiosciencesRRID: AB 39504Mouse monoclonal anti-RBFOX3/NeuNNovusCat# 610774, RRID: AB 398095Mouse monoclonal anti-neuronal-specific enolase (NSE)NovusCat#AB951; RRID: AB 1103747Rabbit polyclonal anti-neuronal-specific enolase (NSE)MilliporeCat#AB951; RRID: AB 22390Rabbit polyclonal anti-TrkB-FL (C-ter region)Santa CruzCat#sc-360 RRID: AB 2227516Rabbit polyclonal anti-TrkB-FL (C-ter region)Santa CruzCat#A-11030; RRID: AB 23254Goat anti-mouse IgG Alexa Fluor 546Molecular ProbesCat#A-11030; RRID: AB 235805Goat anti-mouse IgG Alexa Fluor 647Molecular ProbesCat# A11034 RRID: AB 253540Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11034 RRID: AB 2536217Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11034 RRID: AB 2536217Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A110351 RRID: AB 2576217Goat anti-rabbit IgG-heavy and light chain-HRPBethylCat# A1008P; RRID: AB 143051Donkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A100178; RRID: AB 1211460Chemical and biological productsThermo FisherCat#17504044Ara C (used at 10 µM)Sigma-AldrichCat#01768; Cat# PHG0311DAPI (used at 0.5 or 5 µg/ml)Molecular ProbesCat# PHG0311Human epidermal provt factor recombinant prot. (EGF)Thermo FisherCat#1610-0718Laminin (used at 4 µg/ml)Sigma-Aldrich	Rabbit polyclonal anti- Iba1	Wako	Cat#019-19741
Mouse monoclonal anti-MEF2D BD Biosciences Cat# 610774, RRID:AB 398095 Mouse monoclonal anti-RBFOX3/NeuN Novus Cat#AB951; RRID:AB 1103747 Rabbit polyclonal anti-neuronal-specific enolase (NSE) Millipore Cat#AB951; RRID:AB 92390 Rabbit polyclonal anti N-ter region of RhoGD1 Santa Cruz Cat#se-360 RRID:AB 2227516 Rabbit polyclonal anti-TrkB-FL (C-ter region) Santa Cruz Cat#se-11; RRID:AB 632554 Rabbit polyclonal anti-TrkB-TI (isoform-specific C-ter) Custom made Cat#A-11030; RRID:AB 632554 Goat anti-mouse IgG Alexa Fluor 546 Molecular Probes Cat#A-21236 RRID:AB 2535805 Goat anti-mouse IgG Alexa Fluor 647 Molecular Probes Cat# A-11030; RRID:AB 2535805 Goat anti-rabbit IgG Alexa Fluor 546 Molecular Probes Cat# A11034 RRID:AB 2535805 Goat anti-rabbit IgG Alexa Fluor 546 Molecular Probes Cat# A10034 RRID:AB 2576217 Goat anti-rabbit IgG Alexa Fluor 546 Molecular Probes Cat# A10034 RRID:AB 143051 Donkey anti-rabbit IgG-heavy and light chain-HRP Bethyl Cat#A10078 RRID:AB 10892625 Donkey anti-rabbit IgG-heavy and light chain-HRP Bethyl Cat#A10078; RRID:AB 1211460 Chemical and biological products		() uito	RRID: AB 839504
RRID:AB398095Mouse monoclonal anti-RBFOX3/NeuNNovusRRID:ABRabbit polyclonal anti-neuronal-specific enolase (NSE)MilliporeCat#AB951; RRID:ABRabbit polyclonal anti N-ter region of RhoGDISanta CruzCat#sc-360 RAtbs: ABRabbit polyclonal anti-TrkB-FL (C-ter region)Santa CruzCat#sc-360 RRID:ABRabbit polyclonal anti-TrkB-FL (C-ter region)Santa CruzCat#sc-360 RRID:ABRabbit polyclonal anti-TrkB-TI (isoform-specific C-ter)Custom madeGoat anti-mouse IgG Alexa Fluor 546Molecular ProbesCat#A-11030; RRID:ABGoat anti-mouse IgG Alexa Fluor 546Molecular ProbesCat#A-21236 RRID:ABGoat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11034 RRID:ABGoat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11034 RRID:ABDonkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A120-108P; RRID:ABDonkey anti-mouse IgG-heavy and light chain-HRPBethylCat#A120-108P; RRID:ABB27 scrum free supplementThermo FisherCat#A90-137P; RRID:ABB27 scrum free supplementThermo FisherCat#17504044Recombinant human/murinc/rat brain-derived neurotrophic factor (BDNF) (used at 100 µM)Molecular ProbesCat#A120-629Donkey at 0.5 or 5 µg/ml)Molecular ProbesCat#13256-029Drumesein Autional Autor factor recombinant prot. (EGF)Thermo FisherCat#A120-6311Human enjedermal growth factor recombinant prot. (EGF)Thermo FisherCat#A161-0718Laminin (used at 1	Mouse monoclonal anti-MEF2D	BD Biosciences	Cat# 610774,
Mouse monoclonal anti-RBFOX3/NeuNNovusCat#NBP1-92693SS; RRID:AB_1103747Rabbit polyclonal anti-neuronal-specific enolase (NSE)MilliporeCat#AB951; RRID:AB_92390Rabbit polyclonal anti N-ter region of RhoGDISanta CruzCat#se-360 RRID: AB_2227516Rabbit polyclonal anti-TrkB-FL (C-ter region)Santa CruzCat#se-360 RRID: AB_2227516Rabbit polyclonal anti-TrkB-TI (isoform-specific C-ter)Custom madeGoat anti-mouse IgG Alexa Fluor 546Molecular ProbesCat#A-11030; RRID:AB_235805Goat anti-mouse IgG Alexa Fluor 647Molecular ProbesCat#A-11034 RRID:AB_2535805Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11034 RRID:AB_253605Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11035 RRID: AB_2576217Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11035 RRID: AB_143051Donkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A120-108P; RRID: AB_1211460Chemical and biological products			RRID:AB_398095
RRID:AB1103747Rabbit polyclonal anti-neuronal-specific enolase (NSE)MilliporeCat#AB951; RRID:ABRabbit polyclonal anti N-ter region of RhoGDISanta CruzCat#sc-360 RRID:ABRabbit polyclonal anti-TrkB-FL (C-ter region)Santa CruzCat#sc-11; RRID:ABRabbit polyclonal anti-TrkB-FL (C-ter region)Santa CruzCat#sc-11; RRID:ABGoat anti-mouse IgG Alexa Fluor 546Molecular ProbesCat#A-11030; RRID:ABGoat anti-mouse IgG Alexa Fluor 546Molecular ProbesCat#A-21236 RRID:ABGoat anti-rabbit IgG Alexa Fluor 488Molecular ProbesCat# A11034 RRID:ABGoat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11035 RRID:ABDonkey anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11035 RRID:ABDonkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A100108P; RRID:ABDonkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A100108P; RRID:ABDonkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A102108P; RRID:ABDonkey anti-mouse IgG-heavy and light chain-HRPBethylCat#A102108P; RRID:ABDonkey anti-mouse IgG-heavy and light chain-HRPCat#A102Donkey anti-mouse IgG-heavy and light chain-HRPCat#A100Donkey anti-mouse IgG-	Mouse monoclonal anti-RBFOX3/NeuN	Novus	Cat#NBP1-92693SS;
Rabbit polyclonal anti-neuronal-specific enolase (NSE)MilliporeCat#AB951; RRID:AB 92390Rabbit polyclonal anti N-ter region of RhoGD1Santa CruzCat#sc-360 RRID:AB 9227516Rabbit polyclonal anti-TrkB-FL (C-ter region)Santa CruzCat#sc-11; RRID:AB 632554Rabbit polyclonal anti-TrkB-T1 (isoform-specific C-ter)Custom madeCat#sc-11; RRID:AB 632554Goat anti-mouse IgG Alexa Fluor 546Molecular ProbesCat#A-11030; RRID:AB 144695Goat anti-mouse IgG Alexa Fluor 647Molecular ProbesCat#A-21236 RRID:AB 144695Goat anti-rabbit IgG Alexa Fluor 488Molecular ProbesCat#A11034 RRID:AB 2576217Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11035 RRID: AB 2576217Goat anti-rabbit IgG-heavy and light chain-HRPBethylCat#A1035 RRID:AB 148051Donkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A10037P; RRID:AB 10892625Donkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A1037P; RRID:AB 10892625Donkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A10768; RTID:AB 1211460Chemical and biological productsThermo FisherCat#17504044Recombinant human/murinc/rat brain-derived neurotrophic factor (BDNF) (used at 100 mg/ml)Molecular ProbesCat#21326-029DAPI (used at 0.5 or 5 µg/ml)Molecular ProbesCat#413256-029Human epidermal growth factor recomb. prot. (bFGF)Thermo FisherCat#13256-029Fluorescein Avidin DVector LaboratoricsCat#A10-0718Laminin (used at 0.5 mg/ml) <td></td> <td></td> <td>RRID:AB_1103747</td>			RRID:AB_1103747
Rabbit polyclonal anti N-ter region of RhoGDI Santa Cruz Cat#sc-360 RRID: AB 2227516 Rabbit polyclonal anti-TrkB-FL (C-ter region) Santa Cruz Cat#sc-310 RRID: AB 2227516 Rabbit polyclonal anti-TrkB-FL (C-ter region) Santa Cruz Cat#sc-11; RRID:AB 25254 Rabbit polyclonal anti-TrkB-FL (isoform-specific C-ter) Custom made RRID:AB 144695 Goat anti-mouse IgG Alexa Fluor 546 Molecular Probes Cat#A-11030; RRID:AB 144695 Goat anti-mouse IgG Alexa Fluor 647 Molecular Probes Cat#A-21236 RRID:AB 2535805 Goat anti-rabbit IgG Alexa Fluor 647 Molecular Probes Cat#A-11034 RRID:AB 2535805 Goat anti-rabbit IgG Alexa Fluor 546 Molecular Probes Cat#A1034 RRID: AB 2535805 Donkey anti-rabbit IgG Alexa Fluor 546 Molecular Probes Cat#A1035 RRID: AB 143051 Donkey anti-rabbit IgG-heavy and light chain-HRP Bethyl Cat#A102-108P; RRID: AB 143051 Donkey anti-mouse IgG-heavy and light chain-HRP Bethyl Cat#A90-137P; RRID:AB 1211460 Chemical and biological products Ara C (used at 10 µM) Sigma-Aldrich Cat#17504044 Recombinant human/murinc/rat brain-derived neurotrophic factor (BDNF) (used at 100 ng/ml) Molecular Probes Cat#17504044 DAPI (used at 0.5 or 5 µg/ml) Molecu	Rabbit polyclonal anti-neuronal-specific enolase (NSE)	Millipore	Cat#AB951;
Rabbit polyclonal anti N-ter region of RhoGDISanta CruzCat#sc-360 RRID: AB 2227516Rabbit polyclonal anti-TrkB-FL (C-ter region)Santa CruzCat#sc-11; RRID:AB 632554Rabbit polyclonal anti-TrkB-T1 (isoform-specific C-ter)Custom madeGoat anti-mouse IgG Alexa Fluor 546Molecular ProbesCat#A-11030; RRID:AB 144695Goat anti-mouse IgG Alexa Fluor 647Molecular ProbesCat#A-21236 RRID:AB 144695Goat anti-rabbit IgG Alexa Fluor 488Molecular ProbesCat# A11034 RRID: AB 2535805Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11034 RRID: AB 2576217Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11035 RRID: AB 12051Donkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A120-108P; RRID: AB 143051Donkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A1037P; RRID: AB 1211460Chemical and biological products			RRID:AB_92390
RRD: AB2227516Rabbit polyclonal anti-TrkB-FL (C-ter region)Santa CruzCat#sc-11; RRID:AB_632554Rabbit polyclonal anti-TrkB-T1 (isoform-specific C-ter)Custom madeRID:AB_632554Goat anti-mouse IgG Alexa Fluor 546Molecular ProbesCat#A-11030; RRID:AB_144695Goat anti-mouse IgG Alexa Fluor 647Molecular ProbesCat#A-21236 RRID:AB_2535805Goat anti-rabbit IgG Alexa Fluor 488Molecular ProbesCat#A11034 RRID:AB_2535805Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11035 RRID:AB_2576217Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11035 RRID:AB_2576217Donkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A120-108P; RRID:AB_10892625Donkey anti-mouse IgG-heavy and light chain-HRPBethylCat#A100-137P; RRID:AB_1211460Chemical and biological productsThermo FisherCat#A17504044Recombinant human/murine/rat brain-derived neurotrophic factor (BDNF) (used at 100 ng/ml)Molecular ProbesCat#D1306DAPI (used at 0.5 or 5 µg/ml)Molecular ProbesCat#D1306Human epidermal growth factor recombinant prot. (EGF)Thermo FisherCat#D1306Human basic fibroblast growth factor recomb. prot. (bFGF)Thermo FisherCat#13256-029Fluorescein Avidin DVector LaboratoricsCat#161-0718Camini (used at 0.5 mg/ml)Sigma-AldrichCat#161-0718Laminin (used at 0.5 mg/ml)Sigma-AldrichCat#12020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#161-0718 <td>Rabbit polyclonal anti N-ter region of RhoGDI</td> <td>Santa Cruz</td> <td>Cat#sc-360</td>	Rabbit polyclonal anti N-ter region of RhoGDI	Santa Cruz	Cat#sc-360
Rabbit polycional anti-TrkB-FL (C-ter region)Santa CruzCat#sc-11; RRID:AB 632554Rabbit polycional anti-TrkB-T1 (isoform-specific C-ter)Custom madeRID:AB 632554Goat anti-mouse IgG Alexa Fluor 546Molecular ProbesCat#A-11030; RRID:AB 144695Goat anti-mouse IgG Alexa Fluor 647Molecular ProbesCat#A-21236 RRID:AB 235805Goat anti-rabbit IgG Alexa Fluor 488Molecular ProbesCat# A11034 RRID:AB 2576217Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11035 RRID: AB 2576217Goat anti-rabbit IgG-heavy and light chain-HRPBethylCat# A10351 RRID: AB 143051Donkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A00-137P; RRID: AB 143051Donkey anti-mouse IgG-heavy and light chain-HRPBethylCat#A00-137P; RRID: AB 1211460Chemical and biological productsAra C (used at 10 µM)Sigma-AldrichCat#C1768; Cat#17504044B27 serum free supplementThermo FisherCat#17504044Recombinant human/murine/rat brain-derived neurotrophic factor (BDNF) (used at 100 ng/ml)Molecular ProbesCat#D1306DAPI (used at 0.5 or 5 µg/ml)Molecular ProbesCat#13256-029Fluorescein Avidin DVector LaboratoriesCat#A 2001Glycine (used at 10 µM)Bio-RadCat#161-0718Laminin (used at 4 µg/ml)Sigma-AldrichCat#12200MTT (used at 0.5 or 5 µg/ml)Molecular ProbesCat#A 2001Guest A 100Sigma-AldrichCat#A2001Guest A 100Sigma-AldrichCat#A2001Guest A 1			RRID: AB_2227516
Rabbit polyclonal anti-TrkB-T1 (isoform-specific C-ter)Custom madeGoat anti-mouse IgG Alexa Fluor 546Molecular ProbesCat#A-11030; RRID:AB 144695Goat anti-mouse IgG Alexa Fluor 647Molecular ProbesCat#A-21236 RRID:AB 2535805Goat anti-rabbit IgG Alexa Fluor 488Molecular ProbesCat# A11034 RRID: AB 2576217Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11034 RRID: AB 2576217Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11035 RRID: AB 143051Donkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A120-108P; RRID: AB 143052Donkey anti-mouse IgG-heavy and light chain-HRPBethylCat#A90-137P; RRID: AB 1211460Chemical and biological productsThermo FisherCat#17504044Ara C (used at 10 µM)Sigma-AldrichCat#21768; B27 serum free supplementB27 serum free supplementThermo FisherCat#17504044Recombinant human/murine/rat brain-derived neurotrophic factor (BDNF) (used at 100 ng/ml)Molecular ProbesCat#13256-022DAPI (used at 0.5 or 5 µg/ml)Molecular ProbesCat#13256-029Fluorescein Avidin DHuman epidermal growth factor recomb. prot. (bFGF)Thermo FisherCat#13256-029Fluorescein Avidin DGlycine (used at 10 µM)Bio-RadCat#161-0718Cat#12020Luminin (used at 4 µg/ml)Sigma-AldrichCat#12020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#A2001Guyene (used at 10 µM)Sigma-AldrichCat#A2001	Rabbit polyclonal anti-TrkB-FL (C-ter region)	Santa Cruz	Cat $\#$ sc-11;
Rabbit polycional and - FrkB-Fr (isofold)-specific C-ter)Cusion madeGoat anti-mouse IgG Alexa Fluor 546Molecular ProbesCat#A-11030; RRID:AB 144695Goat anti-mouse IgG Alexa Fluor 647Molecular ProbesCat#A-21236 RRID:AB 2535805Goat anti-rabbit IgG Alexa Fluor 488Molecular ProbesCat#A 11034 RRID:AB 2576217Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11035 RRID: AB 2576217Goat anti-rabbit IgG-heavy and light chain-HRPBethylCat#A120-108P; RRID:AB 143051Donkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A90-137P; RRID:AB 1211460Chemical and biological products	Dabbit naturalanal anti Tuly D T1 (icafarm anasifia (tan)	Custom modo	KKID:AB_032354
Goat anti-mouse IgG Alexa Fluor 546Molecular ProbesCat#A-11030; RRID:AB_144695Goat anti-mouse IgG Alexa Fluor 647Molecular ProbesCat#A-21236 RRID:AB_2535805Goat anti-rabbit IgG Alexa Fluor 488Molecular ProbesCat# A11034 RRID: AB_2576217Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11035 RRID: AB_2576217Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11035 RRID: AB_143051Donkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A120-108P; RRID: AB_143051Donkey anti-mouse IgG-heavy and light chain-HRPBethylCat#A90-137P; RRID: AB_1211460Chemical and biological productsThermo FisherCat#17504044Recombinant human/murine/rat brain-derived neurotrophic factor (BDNF) (used at 100 ng/ml)PeproTechCat#450-02DAPI (used at 0.5 or 5 µg/ml)Molecular ProbesCat#13256-029Cat#13256-029Human epidermal growth factor recomb. prot. (EGF)Thermo FisherCat#13256-029Fluorescein Avidin DVector LaboratoriesCat#A 2001Glycine (used at 10 µM)Bio-RadCat#161-0718Laminin (used at 4 µg/ml)Sigma-AldrichCat#A 2001Glycine (used at 10 µM)Sigma-AldrichCat#A 2001Glycine (used at 0.5 mg/ml)Sigma-AldrichCat#A 2001Glycine (used at 10 µM)Bio-RadCat#161-0718Laminin (used at 4 µg/ml)Sigma-AldrichCat#A2000MTT (used at 0.5 mg/ml)Sigma-AldrichCat#A2000	Rabbit porycional anti- IIRD-II (isoform-specific C-ter)		
RRID:AB 144695Goat anti-mouse IgG Alexa Fluor 647Molecular ProbesCat#A-21236 RRID:AB 2535805Goat anti-rabbit IgG Alexa Fluor 488Molecular ProbesCat# A11034 RRID: AB 2576217Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11035 RRID: AB 143051Donkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A120-108P; RRID: AB 10892625Donkey anti-mouse IgG-heavy and light chain-HRPBethylCat#A120-108P; RRID: AB 10892625Donkey anti-mouse IgG-heavy and light chain-HRPBethylCat#A120-137P; RRID:AB 10892625Donkey anti-mouse IgG-heavy and light chain-HRPBethylCat#A120-137P; RRID:AB 10892625Donkey anti-mouse IgG-heavy and light chain-HRPBethylCat#A120-108P; RRID:AB 1211460Chemical and biological productsFormo FisherCat#17504044Recombinant human/murine/rat brain-derived neurotrophic factor (BDNF) (used at 100 ng/ml)PeproTechCat#450-02DAPI (used at 0.5 or 5 µg/ml)Molecular ProbesCat#D1306Human epidermal growth factor recombinant prot. (EGF)Thermo FisherCat#13256-029Fluorescein Avidin DVector LaboratoriesCat#A 2001Glycine (used at 10 µM)Bio-RadCat#161-0718Laminin (used at 4 µg/ml)Sigma-AldrichCat#42000MTT (used at 0.5 mg/ml)Sigma-AldrichCat#42000Glycine (used at 10 µM)Sigma-AldrichCat#42000Glycine (used at 10 µM)Sigma-AldrichCat#42000Glycine (used at 10 µM)Sigma-AldrichCat#42020Gly	Goat anti-mouse IgG Alexa Fluor 546	Molecular Probes	Cat#A-11030:
Goat anti-mouse IgG Alexa Fluor 647Molecular ProbesCat#A-21236 RRID:AB 2535805Goat anti-rabbit IgG Alexa Fluor 488Molecular ProbesCat# A11034 RRID: AB 2576217Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11035 RRID: AB 143051Donkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A120-108P; RRID:AB 10892625Donkey anti-mouse IgG-heavy and light chain-HRPBethylCat#A120-108P; RRID:AB 1211460Chemical and biological productsThermo FisherCat#17504044Recombinant human/murine/rat brain-derived neurotrophic factor (BDNF) (used at 100 ng/ml)PeproTechCat#17504044DAPI (used at 100 ng/ml)Molecular ProbesCat#11504044DAPI (used at 100 ng/ml)Molecular ProbesCat#11504044DAPI (used at 100 ng/ml)Molecular ProbesCat#115060Human epidermal growth factor recombinant prot. (EGF)Thermo FisherCat#13256-029Fluorescein Avidin DVector LaboratoriesCat#A 2001Glycine (used at 10 µM)Bio-RadCat#161-0718Laminin (used at 4 µg/ml)Sigma-AldrichCat#2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#2020			RRID:AB 144695
RRID: AB2535805Goat anti-rabbit IgG Alexa Fluor 488Molecular ProbesCat# A11034 RRID: ABGoat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11035 RRID: ABDonkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A120-108P; RRID: ABDonkey anti-mouse IgG-heavy and light chain-HRPBethylCat#Cat#A120-108P; RRID: ABDonkey anti-mouse IgG-heavy and light chain-HRPBethylCat#A120-108P; RRID: ABDonkey anti-mouse IgG-heavy and light chain-HRPSigma-AldrichCat#Cat#Co2Cat#a C (used at 10 µM)Bio-RadCat#H161-0718Laminin (used at 4 µg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5	Goat anti-mouse IgG Alexa Fluor 647	Molecular Probes	Cat#A-21236
Goat anti-rabbit IgG Alexa Fluor 488Molecular ProbesCat# A11034 RRID: AB 2576217Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11035 RRID: AB 143051Donkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A120-108P; RRID: AB 10892625Donkey anti-mouse IgG-heavy and light chain-HRPBethylCat#A90-137P; RRID: AB 1211460Chemical and biological productsKara C (used at 10 μM)Sigma-AldrichCat#C1768; Cat#17504044B27 serum free supplementThermo FisherCat#450-02 Cat#450-02Cat#450-02factor (BDNF) (used at 100 ng/ml)Molecular ProbesCat#D1306DAPI (used at 0.5 or 5 μg/ml)Molecular ProbesCat#A123256-029Fluorescein Avidin DVector LaboratoriesCat#A2001Glycine (used at 10 μM)Bio-RadCat#A10306Human basic fibroblast growth factor recomb. prot. (bFGF)Thermo FisherCat#A2001Glycine (used at 10 μM)Bio-RadCat#A2001Glycine (used at 10 μM)Sigma-AldrichCat#A2001Glycine (used at 10 μM)Bio-RadCat#A2001Glycine (used at 10 μM)Sigma-AldrichCat#A2001Glycine (used at 10 μM)Sigma-AldrichCat#A2001Glycine (used at 10 μM)Sigma-AldrichCat#A2001Glycine (used at 0.5 mg/ml)Sigma-AldrichCat#A5555Mitter (used at 0.5 mg/ml)Sigma-AldrichCat#A5555			RRID:AB 2535805
RRID: AB_2576217Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11035 RRID: AB 143051Donkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A120-108P; RRID: AB 10892625Donkey anti-mouse IgG-heavy and light chain-HRPBethylCat#A120-137P; RRID: AB_1211460Chemical and biological productsCat#A90-137P; RRID: AB_1211460Ara C (used at 10 µM)Sigma-AldrichCat#C1768; Cat#17504044B27 serum free supplementThermo FisherCat#450-02 factor (BDNF) (used at 100 ng/ml)DAPI (used at 0.5 or 5 µg/ml)Molecular ProbesCat#D1306Human epidermal growth factor recombinant prot. (EGF)Thermo FisherCat#13256-029Fluorescein Avidin DVector LaboratoriesCat#A 2001Glycine (used at 10 µM)Bio-RadCat#A 2001Glycine (used at 10 µM)Sigma-AldrichCat#A 2001Glycine (used at 10 µM)Bio-RadCat#A 2001Glycine (used at 10 µM)Sigma-AldrichCat#A 2001Glycine (used at 10 µM)Bio-RadCat#A 2001Glycine (used at 10 µM)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#L2020	Goat anti-rabbit IgG Alexa Fluor 488	Molecular Probes	Cat# A11034
Goat anti-rabbit IgG Alexa Fluor 546Molecular ProbesCat# A11035 RRID: AB 143051Donkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A120-108P; RRID:AB 10892625Donkey anti-mouse IgG-heavy and light chain-HRPBethylCat#A90-137P; RRID:AB 10892625Donkey anti-mouse IgG-heavy and light chain-HRPBethylCat#C1768; RRID:AB 1211460Chemical and biological productsAra C (used at 10 μM)Sigma-AldrichCat#C1768; Cat#C1768;B27 serum free supplementThermo FisherCat#C1768; Cat#450-02factor (BDNF) (used at 100 ng/ml)Molecular ProbesCat#D1306Human epidermal growth factor recombinant prot. (EGF)Thermo FisherCat#13256-029Fluorescein Avidin DVector LaboratoriesCat#A 2001Glycine (used at 10 μM)Bio-RadCat#161-0718Laminin (used at 4 μg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#L2020			RRID: AB_2576217
RRID: AB 143051Donkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A120-108P; RRID:AB 10892625Donkey anti-mouse IgG-heavy and light chain-HRPBethylCat#A90-137P; RRID:AB_1211460Chemical and biological productsCat#A90-137P; RRID:AB_1211460Ara C (used at 10 µM)Sigma-AldrichCat#C1768;B27 serum free supplementThermo FisherCat#17504044Recombinant human/murine/rat brain-derived neurotrophic factor (BDNF) (used at 100 ng/ml)PeproTechCat#450-02DAPI (used at 0.5 or 5 µg/ml)Molecular ProbesCat#D1306Human epidermal growth factor recombinant prot. (EGF)Thermo FisherCat#12256-029Fluorescein Avidin DVector LaboratoriesCat#A 2001Glycine (used at 10 µM)Bio-RadCat#161-0718Laminin (used at 4 µg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#L2020	Goat anti-rabbit IgG Alexa Fluor 546	Molecular Probes	Cat# A11035
Donkey anti-rabbit IgG-heavy and light chain-HRPBethylCat#A120-108P; RRID:AB_10892625Donkey anti-mouse IgG-heavy and light chain-HRPBethylCat#A90-137P; RRID:AB_1211460Chemical and biological productsAra C (used at 10 μM)Sigma-AldrichCat#C1768;Ara C (used at 10 μM)Sigma-AldrichCat#17504044Recombinant human/murine/rat brain-derived neurotrophic factor (BDNF) (used at 100 ng/ml)PeproTechCat#450-02DAPI (used at 0.5 or 5 μg/ml)Molecular ProbesCat#D1306Human basic fibroblast growth factor recombinant prot. (EGF)Thermo FisherCat#13256-029Fluorescein Avidin DVector LaboratoriesCat#A 2001Glycine (used at 10 μM)Bio-RadCat#161-0718Laminin (used at 4 μg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#L2020		D 4 1	RRID: AB 143051
Donkey anti-mouse IgG-heavy and light chain-HRPBethylCat#A90-137P; RRID:AB_1211460Chemical and biological productsRRID:AB_1211460Ara C (used at 10 μM)Sigma-AldrichCat#C1768;B27 serum free supplementThermo FisherCat#17504044Recombinant human/murine/rat brain-derived neurotrophic factor (BDNF) (used at 100 ng/ml)Molecular ProbesCat#D1306DAPI (used at 0.5 or 5 µg/ml)Molecular ProbesCat#D1306Human epidermal growth factor recombinant prot. (EGF)Thermo FisherCat#13256-029Fluorescein Avidin DVector LaboratoriesCat#A 2001Glycine (used at 10 μM)Bio-RadCat#161-0718Laminin (used at 4 µg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#M5655	Donkey anti-rabbit IgG-heavy and light chain-HRP	Bethyl	Cat#A120-108P;
Donkey and fight chain-fixerBethyrCat#A90-13/F, RRID:AB_1211460Chemical and biological productsAra C (used at 10 μM)Sigma-AldrichCat#C1768;B27 serum free supplementThermo FisherCat#17504044Recombinant human/murine/rat brain-derived neurotrophic factor (BDNF) (used at 100 ng/ml)Molecular ProbesCat#D1306DAPI (used at 0.5 or 5 μg/ml)Molecular ProbesCat#D1306Human epidermal growth factor recombinant prot. (EGF)Thermo FisherCat#PHG0311Human basic fibroblast growth factor recomb. prot. (bFGF)Thermo FisherCat#A 2001Glycine (used at 10 μM)Bio-RadCat#I61-0718Laminin (used at 4 μg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#M5655	Dankay anti mayoo LaC haavy and light shain HDD	Dathul	<u>KKID:AB_10892625</u>
RKRD.RB_1211400Chemical and biological productsAra C (used at 10 μM)Sigma-AldrichCat#C1768;B27 serum free supplementThermo FisherCat#17504044Recombinant human/murine/rat brain-derived neurotrophic factor (BDNF) (used at 100 ng/ml)PeproTechCat#450-02DAPI (used at 0.5 or 5 µg/ml)Molecular ProbesCat#D1306Human epidermal growth factor recombinant prot. (EGF)Thermo FisherCat#13256-029Fluorescein Avidin DVector LaboratoriesCat#A 2001Glycine (used at 10 µM)Bio-RadCat#161-0718Laminin (used at 4 µg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#M5655	Donkey and mouse igo-neavy and light cham-fike	Deulyi	RRID AR 1211/60
Ara C (used at 10 μM)Sigma-AldrichCat#C1768;B27 serum free supplementThermo FisherCat#17504044Recombinant human/murine/rat brain-derived neurotrophic factor (BDNF) (used at 100 ng/ml)PeproTechCat#450-02DAPI (used at 0.5 or 5 μg/ml)Molecular ProbesCat#D1306Human epidermal growth factor recombinant prot. (EGF)Thermo FisherCat#13256-029Fluorescein Avidin DVector LaboratoriesCat#A 2001Glycine (used at 10 μM)Bio-RadCat#161-0718Laminin (used at 4 μg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#M5655	Chemical and biological products		1211700
Ara C (used at 10 μM)Signa-AldrenCat#C1768,B27 serum free supplementThermo FisherCat#17504044Recombinant human/murine/rat brain-derived neurotrophic factor (BDNF) (used at 100 ng/ml)PeproTechCat#450-02DAPI (used at 0.5 or 5 µg/ml)Molecular ProbesCat#D1306Human epidermal growth factor recombinant prot. (EGF)Thermo FisherCat#PHG0311Human basic fibroblast growth factor recomb. prot. (bFGF)Thermo FisherCat#13256-029Fluorescein Avidin DVector LaboratoriesCat#A 2001Glycine (used at 10 μM)Bio-RadCat#161-0718Laminin (used at 4 µg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#M5655	And C (used at 10 uM)	Sigma Aldrich	$C_{a} + \# C_{1769}$
B27 setum free supprementThermo FisherCat#17304044Recombinant human/murine/rat brain-derived neurotrophic factor (BDNF) (used at 100 ng/ml)PeproTechCat#450-02DAPI (used at 0.5 or 5 µg/ml)Molecular ProbesCat#D1306Human epidermal growth factor recombinant prot. (EGF)Thermo FisherCat# PHG0311Human basic fibroblast growth factor recomb. prot. (bFGF)Thermo FisherCat#13256-029Fluorescein Avidin DVector LaboratoriesCat#A 2001Glycine (used at 10 µM)Bio-RadCat#161-0718Laminin (used at 4 µg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#M5655	P27 gamme free supplement	Sigma-Alunch Thormo Fisher	Cat#C1708;
Recombinant numari numerical brain-derived neurotrophicPeproTechCat#430-02factor (BDNF) (used at 100 ng/ml)Molecular ProbesCat#D1306DAPI (used at 0.5 or 5 µg/ml)Molecular ProbesCat# PHG0311Human epidermal growth factor recombinant prot. (EGF)Thermo FisherCat#13256-029Fluorescein Avidin DVector LaboratoriesCat#A 2001Glycine (used at 10 µM)Bio-RadCat#161-0718Laminin (used at 4 µg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#M5655	B2 / setuli free supplement	Denno Tech	Cat#1/304044
DAPI (used at 100 lg/lll)Molecular ProbesCat#D1306Human epidermal growth factor recombinant prot. (EGF)Thermo FisherCat# PHG0311Human basic fibroblast growth factor recomb. prot. (bFGF)Thermo FisherCat#13256-029Fluorescein Avidin DVector LaboratoriesCat#A 2001Glycine (used at 10 μM)Bio-RadCat#161-0718Laminin (used at 4 μg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#M5655	factor (BDNE) (used at 100 mg/ml)	reprotecti	Cal#430-02
DATE (used at 0.5 of 5 μg/ml)Molecular HolesCat#D1500Human epidermal growth factor recombinant prot. (EGF)Thermo FisherCat# PHG0311Human basic fibroblast growth factor recomb. prot. (bFGF)Thermo FisherCat#13256-029Fluorescein Avidin DVector LaboratoriesCat#A 2001Glycine (used at 10 μM)Bio-RadCat#161-0718Laminin (used at 4 μg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#M5655	DAPI (used at $0.5 \text{ or } 5 \text{ µg/ml}$)	Molecular Probes	Cat#D1306
Human opticitial growth factor recombinant plot. (EGF)Internito FisherCat#13256-029Human basic fibroblast growth factor recomb. prot. (bFGF)Thermo FisherCat#13256-029Fluorescein Avidin DVector LaboratoriesCat#A 2001Glycine (used at 10 μM)Bio-RadCat#161-0718Laminin (used at 4 μg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#M5655	Human enidermal growth factor recombinant prot (ECE)	Thermo Fisher	Cat#D1500
Fluorescein Avidin DVector LaboratoriesCat#A 2001Glycine (used at 10 μM)Bio-RadCat#161-0718Laminin (used at 4 μg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#M5655	Human basic fibroblast growth factor recomb prot. (EGF)	Thermo Fisher	Cat#1100311
Glycine (used at 10 µM)Bio-RadCat#161-0718Laminin (used at 4 µg/ml)Sigma-AldrichCat#L2020MTT (used at 0.5 mg/ml)Sigma-AldrichCat#M5655	Fluerescoin Avidin D	Vector Laboratorios	Cat#15250-029
Laminin (used at 4 µg/ml)Bio-RadCat#101-0/18MTT (used at 0.5 mg/ml)Sigma-AldrichCat#M5655NITE (used at 0.5 mg/ml)Sigma-AldrichCat#M5655	Clycine (used at 10 µM)	Rio-Rad	Cat#161-0718
MTT (used at 0.5 mg/ml) Sigma-Aldrich Cat#L2020 NITE (used at 0.5 mg/ml) Sigma-Aldrich Cat#M5655	Laminin (used at $4 \mu a/m$)	Sigma-Aldrich	Cat#I 2020
	MTT (used at 0.5 mg/ml)	Sigma-Aldrich	Cat#L2020
NMDA (used at 100 µM)	NMDA (used at 100 µM)	Tocris	Cat#0114

PhosSTOP phosphatases inhibitor cocktail tablets	Roche	Cat#04 906 837 001
Poly-L-Lysine (used at 100 µg/ml)	Sigma-Aldrich	Cat#P1524
Prolong Diamond antifade reagent	Molecular Probes	Cat#P36970
Complete protease inhibitor cocktail tablets	Roche	Cat#11 697 498 001
Rose Bengal (used at 20 mg/kg)	Sigma-Aldrich	Cat#R3877
Strentavidin resin	GenScript	Cat#L00353
TTC (used at 2%)	Sigma-Aldrich	Cat#T8877
Pontidos		
Bio_sTT1 _{c:} (Biotin_VGRKKPROPRREVIEHKIPIDG)	GenScript	N/A
Dio-ST Tr _{Ct} (Diotin-TONKKINQAINI VELTIKITEDO)	GenSerint	N/A N/A
Bio TT1 _c (Biotin VGRKKPROPRPDEVI EHKIDI DG)	GenScript	N/A N/A
$\mathbf{T}\mathbf{M}_{\text{res}} \left(V C \mathbf{D} K K \mathbf{D} \mathbf{D} \mathbf{D} \mathbf{D} \mathbf{D} \mathbf{D} \mathbf{D} \mathbf{D}$	Conformat	
	GenScript	N/A
TTI _{Ct} (YGRKKRRQRRRPPFVLFHKIPLDG)	GenScript	N/A
Primer sequences for real-time PCR		
GluN1-F (all isoforms): TCCACCAAGAGCCCTTCGTG	Sigma-Aldrich	N/A
GluN1-R (all isoforms): GCCCCGTACAGATCACCTTC	Sigma-Aldrich	N/A
GluN2A-F: ACGACTGGGACTACAGCCTG	Sigma-Aldrich	N/A
GluN2A-R: CTTCTCTGCCTGCCCATAGC	Sigma-Aldrich	N/A
TrkB-FL-F (isof. spec.): TATCTTCACCCACCTCAAAC	Sigma-Aldrich	N/A
TrkB-FL-R (isof. spec.): GAGAGACTTGACCTGAGCAC	Sigma-Aldrich	N/A
TrkB-T1-F (isof. spec.): GGGGCTGTGCTGCTTGGT	Sigma-Aldrich	N/A
TrkB-T1-R (isof. spec.): GCTGCGGACATCTTTGGAGA	Sigma-Aldrich	N/A
BDNF-F (all isoforms): TGGCTGACACTTTTGAGCAC	Sigma-Aldrich	N/A
BDNF-R (all isoforms): GTTTGCGGCATCCAGGTAAT	Sigma-Aldrich	N/A
GAPDH-F: TTGCCATCAACGACCCCTTC	Sigma-Aldrich	N/A
GAPDH-R: GCCTTGACTGTGCCGTTGAA	Sigma-Aldrich	N/A
NSE-F: ACAGAATGGGGCTGTGTACC	Sigma-Aldrich	N/A
NSE-R: TGGCAACTGTGGGACATGGC	Sigma-Aldrich	N/A
Critical Commercial Assays		
BCA Protein Assay Kit	Thermo Fisher	Cat# 23225
Clarity Western ECL Blotting Substrate	BioRad	Cat# 1705060
L inofectamine 2000	Life Technologies	Cat#11668019
Emperimental Madala: Oneganiama/Straina	Life reelihologies	Cut/1100001)
		NT/A
Balb/c inbred mice (Balb/cOlaHsd)	Envigo Rms, Spain S.L.	N/A
wistar Rat embryos (E18)	In site facility	N/A
Recombinant DNA		
pCRE (25-mer oligonucleotide with the sequence of two TrkB CREs subcloned into pTK-Luc)	(Deogracias et al.,2004)	N/A
pMEF2 (pRSRF; -307 to -242 of Nur77 promoter with	(Woronicz et al., 1995)	N/A
two MEF2 sites subcloned into pGL2-basic)		
pMEF2mut (two inactivating point mutations in	(Woronicz et al., 1995)	N/A
each MEF2 site of pMEF2)		
Software		
Adobe InDesign		RRID:SCR_021799
Adobe Photoshop		RRID:SCR_014199
G*Power 3.1.9.7	http://www.gpower.hhu.de/	RRID:SCR_013726
GraphPad Prism		RRID:SCR_002798
ImageJ	https://imagej.net/	RRID:SCR_003070
R Project for Statistical Computing	http://www.r-project.org	RRID:SCR_001905
String	http://string.embl.de/	RRID:SCR_005223
Other		

Tissue-Tek O.C.T Compound	Sakura	Cat#4583
Protran Western blotting nitrocellulose membrane	GE Healthcare	Cat#GE10600002

Table S2. Selected Bio-sTT1_{Ct}- interacting proteins showing significantly altered levels after different treatment conditions

Uniprot ID	Gene name	logFC	p-value	Description	Confidence	n° of unique peptides
A0A8I5Y1J2		2.526	0.012	Peptidyl-prolyl cis-trans isomerase OS=Rattus norvegicus OX=10116 GN=ENSRNOG00000067128 PE-3 SV-1	High	2
A0A8I5Y9S3	Mrpl48	-1.393	0.025	Mitochondrial ribosomal protein L48 OS=Rattus norvegicus OX=10116 GN=Mrpl48 PE=1 SV=1	High	3
A0A8I5ZLV1	Mtx2	2.616	0.031	Metaxin 2 OS=Rattus norvegicus OX=10116 GN=Mtx2 PE=1 SV=1	High	2
A0A8I5ZS37	Gak	1.724	0.028	Cyclin G associated kinase OS=Rattus norvegicus OX=10116 GN=Gak PE=1 SV=1 CDC23 (Cell division cycle 23	High	3
A0A8I5ZTG0	Cdc23	-0.776	0.040	yeast, homolog), isoform CRA_b OS=Rattus norvegicus OX=10116 GN=Cdc23 PE=1 SV=1	High	6
A0A8I6A326	Rbm6	-1.285	0.011	RNA binding motif protein 6 OS=Rattus norvegicus OX=10116 GN=Rbm6 PE=1 SV-1	High	3
A0A8I6A9H6	Rnps1	1.442	0.044	RNA-binding protein with serine- rich domain 1 OS=Rattus norvegicus OX=10116 GN=Rnps1 PE=3 SV=1	High	3
A0A8I6AKZ9	Actr1b	1.323	0.006	Actin related protein 1B OS=Rattus norvegicus OX=10116 GN=Actr1b PE=1 SV=1	High	2
A0A8I6AN44	Nfasc	-1.092	0.027	Neurofascin OS=Rattus norvegicus OX=10116 GN=Nfasc PE=3 SV=1	High	9
A0A8I6GIS7	Relch	-1.207	0.031	RAB11 binding and LisH domain, coiled-coil and HEAT repeat containing OS=Rattus norvegicus OX=10116 GN=Relch PE=1 SV=1	High	2
A0A8L2RBV8	Rp139-ps13	0.864	0.038	Ribosomal protein L39 OS=Rattus norvegicus OX=10116 GN=Rpl39 PE=1 SV=1	High	2
F7EN52	Cyp46a1	2.113	0.027	Cytochrome P450, family 46, subfamily a, polypeptide 1 OS=Rattus norvegicus OX=10116 GN=Cyp46a1 PE=3 SV=3	High	3
P62893	Rpl39	0.864	0.038	60S ribosomal protein L39 OS=Rattus norvegicus OX=10116 GN=RpI39 PE=1 SV-2	High	2
Q5RJT2	Ftsj3	1.666	0.040	pre-rRNA 2'-O-ribose RNA methyltransferase FTSJ3 OS=Rattus norvegicus OX=10116 GN=Ftsj3 PE=1 SV=1	High	3
Q66H79	Trim32	1.694	0.044	Tripartite motif protein 32 OS=Rattus norvegicus OX=10116 GN=Trim32 PE=1 SV=1	High	2

BDNF vs basal conditions

NMDA vs BDNH	F conditions	_				
Uniprot ID	Gene name	logFC	p-value	Description	Confidence	n° of unique peptides
A0A0G2JYD8	Capn5	1.827	0.008	Calpain 5 OS=Rattus norvegicus OX=10116 GN=Capn5 PE=1 SV=2	High	3
A0A8I5ZQN4	Myo1c	0.906	0.035	Myosin 1C OS=Rattus norvegicus OX=10116 GN=Myo1c PE=1 SV=1	High	7
A0A8I6G6A0	Gnall	1.596	0.039	G protein subunit alpha 11 OS=Rattus norvegicus OX=10116 GN=Gna11 PE=1 SV=1	High	4
A0A8L2R9X1	Abhd6	2.091	0.024	Abhydrolase domain containing 6, acylglycerol lipase OS=Rattus norvegicus OX=10116 GN=Abhd6 PE=1 SV=1	High	2
A2VCW8	Septin7	0.863	0.042	Septin OS=Rattus norvegicus OX=10116 GN=Septin7 PE=1 SV=1	High	10
B2GUZ3	Mthfd11	1.816	0.045	formatetetrahydrofolate ligase OS=Rattus norvegicus OX=10116 GN=Mthfd11 PE=1 SV=1	High	5
D4A280	Pak5	2.383	0.046	Serine/threonine-protein kinase PAK 5 OS=Rattus norvegicus OX=10116 GN=Pak5 PE=1 SV=1	High	2

NMDA vs basal conditions

Uniprot ID	Gene name	logFC	p-value	Description	Confidence	n° of unique peptides
A0A1W2Q6D9	Otud6b	1.158	0.027	ubiquitinyl hydrolase 1 OS=Rattus norvegicus OX=10116 GN=Otud6b PE=1 SV=2	High	2
A0A8I5ZN55	Ago2	-1.303	0.048	Protein argonaute-2 OS=Rattus norvegicus OX=10116 GN=Ago2 PE=1 SV=1	High	6
A0A8I5ZNV8	Pcnp	-1.018	0.034	PEST proteolytic signal- containing nuclear protein OS=Rattus norvegicus OX=10116 GN=Pcnp PE=1 SV=1	High	4
A0A8I6A326	Rbm6	-1.589	0.034	RNA binding motif protein 6 OS=Rattus norvegicus OX=10116 GN=Rbm6 PE=1 SV=1	High	3
A0A8J8YFZ4	Ptbp1	-2.049	0.013	Polypyrimidine tract-binding protein 1 OS=Rattus norvegicus OX=10116 GN=Ptbp1 PE=1 SV=1	High	6
A0A8L2R9X1	Abhd6	1.798	0.033	Abhydrolase domain containing 6, acylglycerol lipase OS=Rattus norvegicus OX=10116 GN=Abhd6 PE=1 SV=1	High	2
D4A280	Pak5	2.460	0.038	Serine/threonine-protein kinase PAK 5 OS=Rattus norvegicus OX=10116 GN=Pak5 PE=1 SV=1	High	2

Uniprot ID and gene name are presented for each protein, together with its logFC obtained after comparison of the two indicated conditions. Individual *P*-values for each protein are also shown. Protein description together with quality parameters of detection (confidence and n° of unique peptides) are also shown.

Table S3. *Arhgdia* (RhoGDI1) and the group of selected Bio-TT1_{Ct}-interacting proteins included in the Reactome function "Signalling by RhoGTPases"

Uniprot ID	Gene name	Description	Confidence	n° of unique peptides
A0A0G2K860	Arhgef12	Rho guanine nucleotide exchange factor 12 OS=Rattus norvegicus OX=10116 GN=Arhgef12 PE=1 SV=2	High	3
A0A8I5Y117	Scfd1	Sec1 family domain containing 1 OS=Rattus norvegicus OX=10116 GN=Scfd1 PE=1 SV=1	High	3
A0A8I5YC36	Ranbp2	RAN binding protein 2 OS=Rattus norvegicus OX=10116 GN=Ranbp2 PE=1 SV=1	High	4
A0A8I6AHE2	Dynclil	Dynein cytoplasmic 1 intermediate chain 1 OS=Rattus norvegicus OX=10116 GN=Dync1i1 PE=3 SV=1	High	5
A0A8I6ALK1	Nup98	Nuclear pore complex protein Nup98- Nup96 OS=Rattus norvegicus OX=10116 GN=Nup98 PE=1 SV=1	High	2
A0A8I6ALV8	Tubala	Tubulin alpha chain OS=Rattus norvegicus OX=10116 GN=Tuba1b PE=3 SV=1	High	5
A0A8I6G6K9	Elmo2	Engulfment and cell motility 2 OS=Rattus norvegicus OX=10116 GN=Elmo2 PE=1 SV=1	High	2
A0A8I6GBC3	Actr2	Actin related protein 2 OS=Rattus norvegicus OX=10116 GN=Actr2 PE=1 SV=1	High	7
A0A8I6GLR6	Prex1	Phosphatidylinositol-3,4,5- trisphosphate-dependent Rac exchange factor 1 OS=Rattus norvegicus OX=10116 GN=Prex1 PE=4 SV=1	High	2
A0A8L2QIZ5	C1qbp	Complement component 1 Q subcomponent-binding protein, mitochondrial OS=Rattus norvegicus OX=10116 GN=C1qbp PE=1 SV=1	High	5
A0A8L2R7U3	Tubb2b	Tubulin beta chain OS=Rattus norvegicus OX=10116 GN=Tubb2b PE=3 SV=1	High	2
B2GV73	Arpc3	Actin-related protein 2/3 complex subunit 3 OS=Rattus norvegicus OX=10116 GN=Arpc3 PE=1 SV=1	High	4
D3ZRB3	Rhof	RCG21806 OS=Rattus norvegicus OX=10116 GN=Rhof PE=1 SV=2 Serine/threonine-protein kinase PAK 5	High	2
D4A280	Pak5	OS=Rattus norvegicus OX=10116 GN=Pak5 PE=1 SV=1	High	2
G3V6S0	Sptbn1	Spectrin beta chain OS=Rattus norvegicus OX=10116 GN=Sptbn1 PE=1 SV=4	High	91
P68035	Actc1	Actin, alpha cardiac muscle 1 OS=Rattus norvegicus OX=10116 GN=Actc1 PE=2 SV=1	High	3
Q5XI73	Arhgdia	Rho GDP-dissociation inhibitor 1 OS=Rattus norvegicus OX=10116 GN=Arhgdia PE=1 SV=1	High	8
Q7TT49	Cdc42bpb	Serine/threonine-protein kinase MRCK beta OS=Rattus norvegicus OX=10116 GN=Cdc42bpb PE=1 SV=1	High	24

Uniprot ID and gene name are presented for each protein together with its description and quality parameters of detection (confidence and n° of unique peptides detected).

Supplementary methods

Antibodies, chemicals and recombinant proteins, cell-penetrating peptides, primers, commercial assays, experimental models, recombinant DNAs and software are described in the Table S1, included as Supplementary material.

Experimental models

Animal procedures were performed following European Union Directive 2010/63/ EU and were approved by the CSIC and Comunidad de Madrid (Ref PROEX 276.6/20) ethics committees. Housing facilities were approved by Comunidad de Madrid (#ES 280790000188) and conform to official regulations. Animals had standard health and immune status and were looked after by professional caretakers daily. Male and female mice were kept in groups of up to five in standard IVC cages while 1–2 pregnant rats occupied standard cages, always containing bedding and nesting materials. Animals were under controlled lighting conditions (12 h light cycles), relative humidity and temperature. Irradiated food and water were provided *ad libitum*. We tried to minimize animal suffering and reduce the number of sacrificed animals.

Mice model of ischemia by microvascular photothrombosis to study in vivo excitotoxicity

This model of permanent focal ischemia mimics embolic or thrombotic occlusion of small arteries, which is frequently found in human stroke, and causes a focal brain damage with histological and MRI correlations to human patterns [1]. The procedure was as previously described [2] with some modifications. Male and female Balb/cOlaHsd mice (23-31 g; 10-16 weeks of age; Envigo Rms, Spain S.L.) were allowed to acclimatize to our facilities for at least 1 week after delivery, before ischemic induction. They were kept with the same cage mates, and daily examined to check their health status. Mice were anesthetized with isoflurane (5% for induction, 2% for maintenance in oxygen; Abbot Laboratories, Madrid, Spain) and then placed in a stereotaxic frame (Narishige Group, Tokyo, Japan). Body temperature was maintained at 36–37°C using a self-regulating heating blanket (Cibertec, Madrid, Spain). A midline scalp incision was made, the skull was exposed, and both Bregma and Lambda points were identified. A cold-light (Schott KL 2500 LCD; Schott Glass, Mainz, Germany) with a fibre optic bundle of 1.5 mm in diameter was centred using a micromanipulator on the right side, at 0.2 mm anterior and 2.0 mm lateral (+0.2 AP, +2 ML) relative to Bregma. Afterward, the photosensitive dye Rose Bengal (7.5 mg/ml, prepared in sterile saline; Sigma-Aldrich) was administered by retro-orbital injection of the venous sinus, for intravenous (i.v.) vascular access, to a body dose of 20 mg/kg. Five minutes later, the brain was illuminated (600 lms, 3,000 K) through the intact skull for 10 min. Brain injury encompasses damage to vascular endothelium, platelet activation, and subsequent microvascular thrombotic occlusion of the irradiated region [3]. The areas underneath this stereotaxic position that result irradiated are the primary motor cortex and the primary somatosensory cortex (hindlimb and forelimb). Once completed surgery, the incision was sutured, and mice were allowed to recover.

For neuroprotection experiments, a single dose (3 nmol/g) of peptides TMyc or $TT1_{Ct}$ (> 95% purity; GenScript), solubilized as 2.5 mM solutions in 0.9% NaCl, we retro-orbitally injected 1 h after damage initiation. In order to improve their plasma stability, these peptides are N- and C-ter modified by, respectively, acetyl and amide groups. Mice were not subjected to other procedures before ischemia and were naive to drug or peptide treatment. Animals were randomly allocated to the experimental groups and the researchers doing the experiments were blind respect to treatment. There was not a previous estimation of sample size and experiments were independently carried out as indicated. Twenty-four hours after damage induction, animals completed the beam walking test and were subsequently sacrificed by cervical dislocation to measure the infarct volume, as indicated below. Before that, brains were sectioned into serial 1-mm-thick coronal slices using a mouse brain matrix (Stoelting, Wood Dale, IL, USA), and sections were completely stained with 2% TTC at room temperature to visualize the cortical infarcts. For immunohistochemistry, animals were deeply anesthetized 5 or 24 h after brain damage and intracardially perfused with cold PBS and 4% paraformaldehyde in PBS as indicated below.

Primary cultures of rat cortical neurons to study in vitro excitotoxicity

The dissected cerebral cortices of Wistar rat embryos (E18), both genders included, were mechanically dissociated in culture medium (Minimum Essential Medium, MEM, Life Technologies) followed by seeding of the cell suspension at a density of 1.1×10^6 cells/ml, prepared in the same medium supplemented with 22.2 mM glucose, 0.1 mM glutamax, 5% foetal bovine serum (FBS), 5% donor horse serum (HS), 100 U/ml penicillin and 100 µg/ml streptomycin. Before seeding, the plates were treated overnight at 37°C with 100 µg/ml poly-L-lysine and 4 µg/ml laminin. The cells were incubated at 37°C in an atmosphere of 5% CO₂ and 95% humidity. Generally, glial growth was inhibited after 7 days *in vitro* (DIVs) by adding 10 µM 5 cytosine β-D-arabinofuranoside (AraC) and experimental treatments took place after 12 DIVs. Excitotoxicity was induced by incubation with NMDA (100 µM) and glycine (10 µM) for the indicated times. As indicated, cultures were preincubated with Tat-derived CPPs before NMDA treatment.

Primary cultures of rat cortical astrocytes

Dissected cerebral cortices obtained from Wistar rat E18 as before were mechanically dissociated in Hank's Balanced Salt Solution (HBSS, Thermo Fisher) containing CaCl₂ and MgCl₂, supplemented with 10 mM Hepes. The resulting cell suspension was then seeded in 75 cm² flasks pretreated with 100 μ g/mL poly-L-lysine at a density of 4x10⁶ cells/ml in proliferation medium. This medium consisted of Dulbecco's Modified Eagle Medium/Nutrient Mixture F-12 with GlutaMAX (DMEM/F-12/GlutaMAX, Thermo Fisher), 10% FBS, 5% HS, 100 U/ml penicillin, 100 μ g/ml streptomycin, 2% B27 supplement, 10 ng/ml epidermal growth factor (EGF), and 10 ng/ml basic fibroblast growth factor (bFGF).

After 4 DIVs, when cultures reached approximately 70–80% confluence, the medium was removed, and the flasks were vigorously washed with PBS to eliminate potential oligodendrocyte contaminants. The cells were then trypsinized using TrypLE Express (Thermo Fisher) and reseeded in flasks under the same conditions at a density of $3x10^6$ cells/ml in proliferation medium. Six days later, astrocytes were washed, trypsinized as before, and seeded at a density of $0.7x10^6$ cells/ml in poly-L-lysine-coated plates for 24 hours. After this period, the plates were washed with PBS, and the medium was replaced with differentiation medium consisting of DMEM/F-12/GlutaMAX, 10% FBS, 100 U/ml penicillin, and 100 µg/ml streptomycin.

Assessment of in vivo and in vitro excitotoxic injury

Measurement of infarct volume

Sections stained with 2% TTC, as before indicated, were subsequently scanned by both rostral and caudal sides. These images were analysed using ImageJ software by an observer blinded to experimental groups. After image calibration, delineated areas of ipsilateral and contralateral hemispheres, and the infarcted region (unstained area) were measured. Considering slices thickness, the corresponding volumes were calculated and corrected for oedema's effect, estimated by comparing total volumes of hemispheres. The corrected infarct volumes were expressed as percentage relative to the contralateral hemisphere, to correct for normal size differences between different animals. For each animal, the mean of results obtained for rostral and caudal sides was calculated.

Beam walking test

Motor coordination and balance were evaluated in mice right before and 24 h after the ischemic insult by measuring the number of contralateral hind paw slips in the beam walk apparatus.

Mice had to walk through a narrow beam $(1 \text{ m} \times 1 \text{ cm} \times 1 \text{ cm})$ placed 50 cm above the tabletop, going from an aversive stimulus (60-W light bulb) to a black goal box with nesting material. Slips taking place in a previously selected central beam segment (50 cm long) were counted. Before damage induction, mice were allowed to cross the beam once, to get acquainted with the test, which they repeated 24 h after photothrombosis, immediately before sacrifice.

Assessment of neuronal injury in cortical cultures

Thiazolyl blue formazan (MTT) reduction assay was used to determine cell viability. At the end of the different treatments, 0.5 mg/ml MTT was added to the medium, and after 2 h of incubation at 37°C the formazan salts formed were solubilized using DMSO. The results were quantified spectrophotometrically at 570 nm. As primary cortical cultures are formed by neurons and glial cells, we determined the magnitude of glia viability on the total values by applying 400 μ M NMDA and 10 μ M glycine to sister cultures 24 h before MTT assay. These conditions induce a complete neuronal death but no glial damage. Thus, after subtracting this absorbance value, the results only represent the viability of neurons present in cultures. Each independent experiments included technical triplicates for every treatment condition. In addition, a minimum of 4 completely independent experiments were performed and analysed as indicated in figure legends. When indicated, the provided viability data correspond to total cell values (neurons and glial cells.

Synthetic peptides

Synthetic cell-penetrating peptides (> 95% purity; GenScript) were used for treatment of primary cortical cultures and mice. All these CPPs contain a HIV-1 Tat sequence (11 aa) linked to a specific TrkB-T1 or c-Myc sequence as indicated (Table S1). Three of these peptides (Bio-TMyc, Bio-sTT1_{Ct} and Bio-TT1_{Ct}) contain a biotin molecule at the N-ter, which serves several purposes: peptide labelling, stabilization, and competence to isolate proteins interacting with TrkB-T1 taking place inside cultured cells. In the other two peptides (TMyc and TT1_{Ct}), the N-ter is acetylated while all five CPPs have their C-ter amidated. The concentration of use is indicated for each experiment.

Pull down assays

Cultures incubated with biotinylated peptides and treated as indicated were lysed at 4°C in NP-40 lysis buffer (1% NP-40, 20 mM Tris HCl pH 8.0, 80 mM NaCl, 20 mM EDTA) containing protease and phosphatase inhibitors (Table S1), conditions that do not disrupt peptide-protein interactions. After establishing protein concentration (BCA Protein Assay Kit), equal protein amounts were combined with streptavidin-agarose beads previously washed with NP-40 lysis buffer and incubated at room temperature for 1.5 h with agitation. Sedimented beads were washed 3 times with the NP-40 lysis buffer to remove possible non-specifically bounded proteins, and peptides and their interacting proteins were released by incubation at 50°C for 40 min in RIPA modified buffer (50 mM Tris HCl pH 8.0, 150 mM NaCl, 1% sodium deoxycholate, 1% NP-40) combined with protease and phosphatase inhibitors as before.

Proteomic analysis of Bio-sTT1_{Ct}-interacting proteins

Proteomic analysis was performed in ProstaR (1.34.6) software tool [4], using the protocol for protein-level analysis [5]. Proteomic abundances from 4 independent experiments were log-transformed and filtered following quality parameters: proteins present in at least >50% of one of the conditions and identified by at least 2 unique peptides. After, data was normalized by LOESS (locally estimated scatterplot smoothing, k=0.7), and POV (partially observed values) and MEC (missing in entire condition values) imputations were made using, respectively, KNN (k-nearest-neighbour k=10) and det.quantile (quantile imputation using quantile = 2,5). Differential analysis was performed using limma (3.16). Multiple comparison p-values where calibrated using slim algorithm and the FDR (false discovery rate) is indicated for each comparison. For PCA (principal component analysis), filtered, normalized and imputed data was used corrected for batch effect using removeBatchEffect {limma} function. First and second principal components are presented. For volcano plot representation ggplot2 (3.5.1) software was used, and for heatmap presentation ComplexHeatmap (2.18.0) software. Enrichment analysis for GO (Gene Ontology terms) and Reactome (Reactome.db) pathway analysis the string.db (12.0) tool was used.

Immunoblot analysis

Cells were lysed in RIPA buffer (50 mM Tris-HCl pH 8.0, 150 mM NaCl, 1% sodium deoxycholate, 1% NP-40, 0.1% SDS, 1 mM DTT), supplemented with protease and phosphatase inhibitors as above, for 30 min at 4°C. Then, cell lysates were sonicated using a Bioruptor apparatus and centrifugated at 4°C for 20 min at 10,000 rpm. The protein concentration in lysates was established with BCA Protein Assay Kit (Thermo Fisher) followed by 5 min denaturalization in SDS-sample buffer at 95°. Equal amounts of total lysates were resolved in Tris-Glycine SDS-PAGE and transferred to a nitrocellulose membrane (GE Healthcare) in 25 mM Tris HCl pH 8.3, 250 mM glycine and 10% methanol, using an electric current of 400 mA for 70 min. After transfer, membranes were stained with a Ponceau S solution to check for its efficacy. Then, membranes were blocked for 30 min with a 5% non-fat

dry milk solution in Tris Buffered Saline-Tween (TBS-T, 20 mM Tris HCl pH 7.5, 137 mM NaCl, 0.05% Tween-20) and incubated overnight at 4°C with primary antibodies. Next, membranes were washed with TBS-T and incubated with the appropriate anti-rabbit or antimouse peroxidase-conjugated secondary antibodies (Bethyl) for 1 h. To conclude, immunoreactivity was detected using Clarity Western ECL Blotting Substrate (BioRad) and band intensity was quantified by densitometric analysis (Photoshop, Adobe). Levels of the protein of interest were normalized using those of neuron-specific enolase (NSE) present in the same sample and expressed relative to values obtained in their respective controls, arbitrarily given a 100% value. NSE was used as a neuronal loading control since it is not affected by NMDA treatment. Exceptionally, levels of TrkB-T1 were normalized to those of β-actin in the comparison of primary mixed cultures and those of astrocytes. Multiple independent experiments were carried out and quantitated as detailed in the figure legends.

Cell transfection and gene reporter assays

Plasmids contained minimal CREB or MEF2 response elements upstream firefly luciferase reporter gene (respectively, pCRE or pMEF2; see details in Supplementary Table 1). Primary cultures, without the AraC treatment, were transfected at 11 DIVs with 0.4 µg of plasmids using Lipofectamine 2000. DNA-liposomes complexes were prepared according to the manufacturer instructions in neurobasal medium (Thermo Fisher Scientific) with 1 mM glutamax (Gibco, Thermo Fisher Scientific) and added to the cell culture. After 2 h of transfection, this mix was replaced with previously collected conditioned medium and incubation proceeded to complete 24 h. Protein extracts were obtained using Passive Lysis Buffer (Promega, Cat# E1941) and luciferase activity was detected by a luminometer (GloMax 96 microplate luminometer, Promega), using REAP buffer (25 mM glycylglycine, 15 mM SO₄Mg, 4 mM EGTA, 15 mM potassium phosphate pH 7.8, 3.3 mM ATP, 1 mM DTT, 75 µM luciferin). Each experiment included technical quadruplicates for every condition, and 5-7 totally independent experiments were analysed as indicated.

RNA extraction and qPCR assay

Total RNA was extracted using QIAcube technology and treated with DNases before cDNA synthesis by the "High Capacity cDNA Reverse Transcription Kit" (Applied Biosystems). A 7900HT Fast real-time PCR system (Applied Biosystems) was used for SYBR green gene expression assays, with the primers indicated in Supplementary Table 1. PCR conditions were 10 min at 95°C, followed by 40 cycles of 15 s at 95°C and 60 s at 60°C. For each independent experiment, we made a specific standard curve for every gene and technical triplicates were

prepared for every sample. Then, we performed a total of 5 completely independent experiments. Results were normalised to the levels of housekeeping genes as indicated in the "Methods" section.

Peptide visualization in primary cultures and immunocytochemistry

Primary cultures, seeded at half the concentration previously indicated, were grown for 13 DIVs on coverslips previously treated overnight at 37°C with poly-L-lysine and laminin as before. Next, cells were incubated with Bio-TMyc or Bio-sTT1_{Ct} (25 μ M, 1 h) or left untreated, and then fixed for 30 min with 4% paraformaldehyde in PBS. The fixed cells were washed several times with PBS to eliminate the remaining paraformaldehyde before they were blocked and permeabilized for 4 h with 4% goat serum and 0.5% Triton X-100 prepared in PBS. Afterwards coverslips were incubated overnight at 4°C with the corresponding primary antibodies, diluted in the blocking solution. After washing with PBS, the coverslips were incubated with the corresponding secondary antibodies together with Fluorescein Avidin D (200 μ g/ml) for 1 h at room temperature and then 10 min with DAPI, to stain DNA. Finally, the samples were mounted in Prolong Diamond. The images were acquired using an inverted Zeiss LSM 710 laser confocal microscope (Jena, Germany) with a 40x or 63x Plan-Apochromatic oil immersion objective. Images correspond to single sections and were normalized for each colour separately and processed for presentation with ImageJ (NIH Image) and Fiji.

Peptide visualization and immunohistochemistry in brain cortex

For peptide visualization in undamaged or damaged brains, animals were deeply anesthetized and intracardially perfused with cold PBS and 4% paraformaldehyde in PBS. Brains were postfixed in the same fixative at 4°C for 24 h and cryoprotected in 30% sucrose for 48 h at 4°C. Coronal frozen sections (30 μ m thick) obtained using a cryostat (Leica, Heidelberg, Germany) were incubated in blocking solution (10% goat serum, 0.5% Triton X-100 in PBS) for 3 h at room temperature, followed by Fluorescein Avidin D (200 μ g/ml), and DAPI (5 μ g/ml) prepared in 4% goat serum for 1 h. When indicated, sections were incubated for 3 h with anti-NeuN in 4% goat serum right after blocking and, after washing, with Alexa Fluor 647conjugated antibodies together with Fluorescein Avidin D and DAPI as before indicated. Finally, sections were mounted and dried on slides, and cover slipped with Prolong Diamond. Image acquisition was performed using an inverted laser confocal microscope as before with a 63x or 40x objective Plan-Apochromatic oil immersion objective. Images were processed as described and correspond to single sections. Background was subtracted using vehicle-injected animals.

Immunohistochemistry

The brains of animals injected one hour after the initiation of the photothrombotic damage with vehicle (saline), TMyc or TT1_{Ct} (3 nmol/g) were processed and cryoprotected as before indicated. The duration of damage before sacrifice was 5 h (male and female) or 24 h (male). Coronal frozen sections were incubated with blocking solution as described in flotation for 1 h at room temperature and then overnight with primary antibodies diluted in 4% goat serum and 0.3% Triton X-100 in PBS at 4°C. After washing, slides were incubated 1 h with Alexa Fluor 546-conjugated goat anti-rabbit or AlexaFluor 647-conjugated goat anti-mouse secondary antibodies to visualize the primary antibodies and DAPI (5 µg/ml). Only for Iba1 staining, blocking solution was 2.5% BSA, 5% goat serum, 0.5% Triton X-100 in PBS, and sections were incubated with primary antibody 24 h at RT followed by washing and Alexa Fluor 546conjugated goat anti-rabbit and DAPI incubation for 1 h. Sections were then washed in PBS, mounted on slides, dried for 15 min at 37°C on a hot plate, or air dried overnight and then cover slipped with DPX or Prolong Diamond. Cell Observed images were obtained with Zeiss microscope Observer.Z1 using a 10x Plan-Apochromatic objective. Confocal images were acquired with a 40x or 63x Plan-Apochromatic oil immersion objective and processed for presentation as above described.

References

1. Pevsner PH, Eichenbaum JW, Miller DC, Pivawer G, Eichenbaum KD, Stern A, et al. A photothrombotic model of small early ischemic infarcts in the rat brain with histologic and MRI correlation. J Pharmacol Toxicol Methods. 2001; 45: 227-33.

2. Tejeda GS, Esteban-Ortega GM, San Antonio E, Vidaurre OG, Díaz-Guerra M. Prevention of excitotoxicity-induced processing of BDNF receptor TrkB-FL leads to stroke neuroprotection. EMBO Mol Med. 2019; 11:e9950.

3. Watson BD, Dietrich WD, Busto R, Wachtel MS, Ginsberg MD. Induction of reproducible brain infarction by photochemically initiated thrombosis. Ann Neurol. 1985; 17: 497-504.

4. Wieczorek S, Combes F, Lazar C, Giai Gianetto Q, Gatto L, Dorffer A, et al. DAPAR & ProStaR: software to perform statistical analyses in quantitative discovery proteomics. Bioinformatics. 2017; 33: 135-6.

5. Wieczorek S, Combes F, Borges H, Burger T. Protein-Level Statistical Analysis of Quantitative Label-Free Proteomics Data with ProStaR. Methods Mol Biol. 2019; 1959: 225-46.