

1 Supplementary Materials for

2 **Blocking TRPM4 alleviates pancreatic acinar cell damage via an NMDA receptor-dependent pathway in**  
3 **acute pancreatitis**

4 Yifan Ren, Qing Cui, Wuming Liu, Hangcheng Liu, Tao Wang, Hongwei Lu, Yi Lv\*, Rongqian Wu\*

5 \*Corresponding author: [rwu001@mail.xjtu.edu.cn](mailto:rwu001@mail.xjtu.edu.cn) or [luyi169@mail.xjtu.edu.cn](mailto:luyi169@mail.xjtu.edu.cn)

6  
7  
8 **This file includes:**

9  
10 Supplementary Materials and Methods

11 Supplementary Table;

12 Figure S1 to Figure S7;

13 Supplementary Figure Legends;

14 Raw western blot images;

15 Report of *trpm4* gene knockout mice.

16 \_\_\_\_\_

17 **Supplementary Materials and Methods**

18 **Trpm4 knockout mouse and genotyping strategy:** Name: C57BL/6NCya-Trpm4 em1/Cya. Serial Number:  
19 KOCMP-68667-Trpm4-B6N-VA. Type: conventional knockout. Primers 1 (Annealing Temperature 60.0 °C), F1:  
20 5'-AGTCCCTGCTCATTACTCTGGG-3', R1: 5'-GAGAAGGTAAAGAAGCCTGTGC-3', Product size: 838 bp. Primers 2  
21 (Annealing Temperature 60.0°C), F2: 5'-TGGGAAGTGGAGGATGAAGTCAG-3', R1: 5'-  
22 GAGAAGGTAAAGAAGCCTGTGC-3', Product size: 643 bp. Homozygotes: one band with 838 bp; Heterozygotes:  
23 two bands with 838 bp and 643 bp; Wildtype allele: one band with 643 bp.



24

25 *Targeting strategy for trpm4 knockout mouse (Enlarged Figure 4A)*

26

27 **Pathological staining:** Pancreas tissue was obtained. The tissue samples were fixed in 10% neutral buffered  
28 formalin at room temperature for 24 - 48 h and then rinsed. Next, the fixed tissues undergo dehydration in  
29 graded ethanol solutions, clearing in xylene, and infiltration with molten paraffin before being embedded in  
30 paraffin blocks. Paraffin - embedded blocks are sectioned to 4 - 5 µm thickness, stretched on a warm water  
31 bath, and air - dried on glass slides. The HE staining process involves deparaffinization, rehydration,  
32 hematoxylin staining for nuclei, differentiation, blueing, eosin staining for cytoplasm and extracellular matrix,  
33 followed by dehydration, clearing, and mounting. The stained sections were digitally scanned with a biopsy  
34 scanner (Pannoramic MIDI, 3DHISTECH, Hungary).

35

36 **Pathological score:** After HE - stained pancreatic tissue sections were prepared as described above, a semi -

37 quantitative histopathological scoring system was employed to evaluate the degree of pancreatic injury. The  
38 scoring was performed independently by two experienced pathologists who were blinded to the experimental  
39 groups. Acinar cell damage: This was evaluated based on the presence of acinar cell necrosis, vacuolization,  
40 and loss of acinar architecture. A score of 0 indicated no damage, 1 for mild damage (less than 25% of acinar  
41 cells affected), 3 for moderate damage (25 - 50% of acinar cells affected), and 5 for severe damage (more than  
42 50% of acinar cells affected). Inflammatory cell infiltration: The amount of inflammatory cell infiltration in the  
43 pancreatic tissue was assessed. A score of 0 represented no infiltration, 1 for mild infiltration (scattered  
44 inflammatory cells), 3 for moderate infiltration (aggregates of inflammatory cells in some areas), and 5 for  
45 severe infiltration (diffuse and extensive inflammatory cell infiltration). Edema: The degree of interstitial  
46 edema was scored. A score of 0 indicated no edema, 1 for mild edema (slight widening of the interstitium), 3  
47 for moderate edema (obvious separation of tissue structures by fluid), and 5 for severe edema (massive  
48 accumulation of fluid with marked tissue distortion).

49

50 **Water content determination:** We assessed the degree of pancreatic edema by evaluating changes in water  
51 content. Briefly, pancreatic tissues were dried in a 72 °C oven for 2 days. The pancreas water content was  
52 calculated as the wet weight/dry weight ratio.

53

54 **Immunohistochemical staining:** 4% formalin-fixed pancreatic tissue was paraffin embedded and sectioned at  
55 4 - 5 µm. After deparaffinization and rehydration, antigen retrieval was done using heat in citrate buffer.  
56 Endogenous peroxidase was blocked with 3% hydrogen peroxide, followed by blocking non - specific binding  
57 with normal serum. Incubation with primary antibody overnight at 4 °C was followed by a biotinylated  
58 secondary antibody and then a streptavidin - peroxidase complex. Chromogenic development with DAB was

59 followed by hematoxylin counterstaining. Sections were dehydrated, cleared in xylene, and mounted.  
60 Immunohistochemical staining was performed and the stained sections were digitally scanned with a biopsy  
61 scanner (Pannoramic MIDI, 3DHISTECH, Hungary). The primary antibodies are TRPM4 (ab106200, Abcam, USA)  
62 and LY6G (ab238132, Abcam, USA).

63

64 **ATP content determination:** Pancreas tissue or AR42J or mouse PACs homogenate was obtained. Preparation  
65 of tissue or cell samples for detection was carried out in accordance with the specifications of the assay kit  
66 (S0026, Beyotime, China). As for the detection process, the ATP standard solution is first diluted to the  
67 appropriate concentration gradient with the ATP detection lysate. Examples include 0.01, 0.03, 0.1, 0.3, 1, 3,  
68 and 10  $\mu\text{M}$ . In subsequent experiments, the concentration range of the standard can be adjusted appropriately  
69 according to the concentration of ATP in the sample. Then, the appropriate amount of ATP test solution is  
70 prepared in the ratio of 100  $\mu\text{L}$  of ATP test solution required for each sample or standard. Take an appropriate  
71 amount of ATP test reagent and dilute it with ATP test reagent diluent at the ratio of 1:4. The diluted ATP  
72 detection reagent is the working liquid for ATP detection for subsequent experiments. To determine the ATP  
73 content of the sample, first add 100  $\mu\text{L}$  of ATP detection working liquid to the detection tube. Leave at room  
74 temperature for 3-5 min, so that all the background ATP is consumed, thereby reducing the background. Then  
75 add 20  $\mu\text{L}$  of sample or standard in the test tube, quickly mix with a micropipette, at least 2 s after the  
76 luminometer or liquid flash meter to determine the RLU value or CPM. Finally, the concentration of ATP in the  
77 sample was calculated according to the standard curve.

78

79 **Transmission electron microscopy (TEM):** Pancreatic tissue was promptly fixed in 2.5% glutaraldehyde in 0.1  
80 M phosphate buffer (pH 7.4) at 4 °C for 2 - 4 h. After rinsing, they were post - fixed with 1% osmium tetroxide.

81 Dehydration was done using a series of ethanol solutions, followed by infiltration with a mixture of absolute  
82 ethanol and epoxy resin, and then embedding in epoxy resin at 60 °C for polymerization. 60-80 nm ultra-thin  
83 sections of pancreatic samples were stained with uranyl acetate and lead citrate. Pancreatic ultrastructure,  
84 such as mitochondria and endoplasmic reticulum, was evaluated using a transmission electron microscope  
85 (HITACHI HT7700, Hitachi, Japan) by a technical microscopist.

86

87 **TUNEL, HSP60, DHE (Dihydroethidium), MitoTracker and Fluo-3 AM staining:** A TUNEL kit (11684795910,  
88 Roche, Switzerland), a HSP60 (Antibody, AG2237, Beyotime, China) kit, two DHE (Dihydroethidium for *in vivo*,  
89 G1045, Servicebio, China or Dihydroethidium for *in vitro*, S0063, Beyotime, China) kit, a MitoTracker staining  
90 kit (Mitochondrial probe, M7512, Thermo Fisher Scientific, Beijing, China) and a Fluo-3 AM staining kit (Calcium  
91 ion fluorescent probe, S1056, Beyotime, China) were used for TUNEL, HSP60, DHE, MitoTracker and Fluo-3  
92 staining *in vivo* and *in vitro* according to the manufacturers' instructions.

93

94 For TUNEL fluorescent staining, pancreatic tissue samples were fixed in 4% paraformaldehyde, dehydrated,  
95 embedded in paraffin, and sectioned at 4 - 5  $\mu$ m. Tissue sections were deparaffinized and rehydrated. All  
96 samples were permeabilized with 0.1% Triton X - 100 in 0.1% sodium citrate. The TUNEL reaction, using a  
97 commercial kit with a mixture of TdT and fluorescein - labeled dUTP, was carried out at 37 °C for 60 min in the  
98 dark, with negative controls lacking TdT. After washing, samples were counterstained with DAPI, washed again,  
99 mounted with antifade medium, and examined under a fluorescence microscope at various magnifications.

100 The percentage of TUNEL - positive apoptotic cells was calculated by counting green - fluorescent (TUNEL -  
101 positive) and blue - fluorescent (total, stained by DAPI) cells in multiple fields, and data were analyzed with  
102 Image J pro software, presented as mean  $\pm$  SEM.

103

104 For immunofluorescence staining, pancreatic tissues obtained from approved sources were fixed in 4%  
105 paraformaldehyde for 24 h, dehydrated, cleared, embedded in paraffin, and sectioned into 4 - 5  $\mu\text{m}$  slices.  
106 Paraffin - embedded tissue sections underwent deparaffinization, rehydration, and antigen retrieval. All  
107 samples were permeabilized with 0.1% Triton X - 100 and blocked with a buffer containing 5% normal serum.  
108 The primary antibody was added and incubated overnight at 4 °C, with negative controls set up. After washing,  
109 the fluorescently - labeled secondary antibody was added and incubated for 1 - 2 h at room temperature.  
110 Following another washing step, nuclei were counterstained with DAPI. Samples were mounted and observed  
111 under a fluorescence microscope. Images were captured at different magnifications, and the intensity and  
112 distribution of the fluorescent signal were analyzed using Image J pro software.

113

114 For DHE, MitoTracker and Fluo-3 AM Staining, fluorescent dyes were added directly to living cells for staining.  
115 In order to prevent cell death, DAPI was not used to label the nucleus of acinic cells *in vitro*. The fluorescent  
116 dye was incubated for 30-60 min and photographed under a confocal laser microscope. Fluorescence intensity  
117 was evaluated using Image J pro software.

118

119 **Flow cytometry (FCM):** For Fluo-3 AM Staining (Calcium ion fluorescent probe, S1056, Beyotime, China),  
120 fluorescent dyes were added directly to living cells for staining. The fluorescent dye was incubated for 30-60  
121 min and the cells were analyzed with a flow cytometry (ACEA Biosciences, Inc.).

122

123 **Biochemical detection:** Serum amylase and lipase were completed in the automatic biochemical analyzer  
124 (Servicebio, Wuhan, CN) according to the biochemical assay kit instructions (C016-1 and A054-1, Nanjing

125 Jiancheng Bioengineering Institute, CN).

126

127 **Enzyme-linked immunosorbent assay (ELISA):** LDH levels in serum and cell supernate were detected by  
128 corresponding ELISA kits (SEB864Mu and SEB864Ra) from Cloud-Clone Corp, CN. Amylase levels in cell  
129 supernate were measured using corresponding ELISA kits (SEB454Ra) from Cloud-Clone Corp, CN. The mouse  
130 IL-6 ELISA kit (SEA079Mu, Cloud-Clone Corp USCN Life Science, Wuhan, China) and tumor necrosis factor- $\alpha$   
131 (TNF- $\alpha$ ) ELISA kit (SEA133Mu, Cloud-Clone Corp USCN Life Science, Wuhan, China) were used for the detection  
132 of the levels of IL-6 and TNF- $\alpha$  according to the manufacturer's instructions. In brief, the testing process is as  
133 follows: (1) Preparation of standards, reagents and samples before the experiment; (2) Add samples (standard  
134 and sample) 100  $\mu$ L and incubate at 37 °C for 1 h; (3) Suck and discard, add 100  $\mu$ L of detection solution A, and  
135 incubate at 37 °C for 1 h; (4) Wash the board three times; (5) Add 100  $\mu$ L detection solution B and incubate at  
136 37 °C for 30 min; (6) Wash the board 5 times; (7) Add 90  $\mu$ L TMB substrate and incubate at 37 °C for 10 - 20  
137 min; (8) Add 50  $\mu$ L of termination liquid and the absorbance (O.D.) was measured at 450 nm wavelength by  
138 enzyme-labeled instrument (Bio-Rad, California, USA), and the sample concentration was calculated.

139

140 **Detection of antioxidant capability:** The total antioxidant capacity of pancreatic tissue or AR42J cells or mouse  
141 PACs was measured by the ferric ion reducing antioxidant power (FRAP) method using a T-AOC Assay Kit (S0116,  
142 Beyotime, China) according to the manufacturer's protocol. Sample preparation: for cell samples, about 1  
143 million cells were collected, placed in 200  $\mu$ L of cold PBS solution, homogenized or ultrasound to fully break  
144 the cells and release the antioxidants in them, centrifuged at 4 °C for about 12000 g for 5 min, and supernatant  
145 was taken for subsequent determination. For the tissue samples, 100  $\mu$ L of cold PBS solution was added to  
146 every 20 mg of tissue, homogenized or ultrasonic to fully break the tissue and release the antioxidants in it,

147 centrifuged at 4 °C for about 12000g for 5 min, and the supernatant was taken for subsequent determination.

148 Determination of total antioxidant capacity: Add 180 µL of FRAP working fluid to each test hole of the 96-well

149 plate. 5 µL of PBS were added into the blank control hole. Add 5 µL of FeSO<sub>4</sub> standard solution of various

150 concentrations into the standard curve detection hole; 5 µL of various samples or 0.15 - 1.5 mM of Trolox were

151 added into the sample test hole as a positive control. A593 was determined after incubation at 37 °C for 3-5

152 min. The total antioxidant capacity of the sample was calculated according to the standard curve.

153

154 **Western blot analysis:** Pancreatic tissues were lysed in cold RIPA (P0013B, Beyotime, CN). The protein

155 concentration was evaluated with the BCA Protein Assay Kit (P0009, Beyotime, Beijing, China). After gel

156 electrophoresis, the protein was transferred to PVDF membrane (FFP70, Beyotime, Beijing, China) and incubation

157 in blocking solution (3% BSA or 5% skimmed milk) at room temperature. Then the membranes were incubated

158 overnight at 4 °C with the primary antibodies. Primary antibodies were diluted in Primary Antibody Dilution Buffer

159 for Western Blot (P0256, Beyotime, Beijing, China). Membranes were washed and then incubated with specific HRP-

160 conjugated secondary antibodies for 1 h at room temperature. Bands were developed using Digital gel image

161 analysis system (Bio-Rad, California, USA) and the gray values of the bands were quantitatively analyzed by

162 Image J software. The antibodies used in this study are listed in the *Supplementary Table* below.


163

164

**Supplementary Table: Antibodies**

<b>Antibody</b>	<b>Item No</b>	<b>Company and location</b>
HSP60 Rabbit Monoclonal Antibody	AG2237	Beyotime Biotechnology, CN
$\beta$ -Actin Mouse Monoclonal Antibody	AF0003	Beyotime Biotechnology, CN
RIP3 Rabbit mAb	15828	Cell Signaling Technology, Beverly, MA, USA
Anti-Cleaved Caspase-3	ab214430	Abcam, Cambridge, MA, USA
Phospho-eIF-2 $\alpha$ (Ser51) Antibody	AF5803	Beyotime Biotechnology, CN
Anti-PGC-1 alpha	ab191838	Abcam, Cambridge, MA, USA
Bax Rabbit mAb	14796	Cell Signaling Technology, Beverly, MA, USA
eIF2 $\alpha$ (D7D3) XP <sup>®</sup> Rabbit mAb	5324	Cell Signaling Technology, Beverly, MA, USA
IRE1 $\alpha$ (14C10) Rabbit mAb	3294	Cell Signaling Technology, Beverly, MA, USA
CHOP Mouse mAb	2895	Cell Signaling Technology, Beverly, MA, USA
Phospho-IRE-1 $\alpha$ (Ser724) Antibody	AF5842	Beyotime Biotechnology, CN
BiP Antibody	3183	Cell Signaling Technology, Beverly, MA, USA
GRP78 Rabbit Polyclonal Antibody	AF0171	Beyotime Biotechnology, CN
PDI (C81H6) Rabbit mAb	3501	Cell Signaling Technology, Beverly, MA, USA
Anti-TRPM4 Antibody	ab106200	Abcam, Cambridge, MA, USA
PINK1 Rabbit Polyclonal Antibody	Af7755	Beyotime Biotechnology, CN
TRPM4 Polyclonal Antibody	PA5-116483	Thermo Fisher Scientific, ShangHai, CN
Goat anti-Mouse IgG antibody	31430	PIONEER Biotechnology, CN
Goat anti-Rabbit IgG antibody	31460	PIONEER Biotechnology, CN
Anti-LY6G antibody	ab238132	Abcam, Cambridge, MA, USA
MFN2 Rabbit Polyclonal Antibody	AF7473	Beyotime Biotechnology, CN


168 **Figure S1. Low dose of 9-phenanthrol had little protective effect on experimental AP.** (A) Representative  
 169 photos of H&E staining of the pancreas (200X); (B) Pancreatic injury scores; (C) Percentages of necrotic areas.  
 170 n = 6, error bars indicate the SEM; n.s: non-significant; AP, acute pancreatitis; LPS, lipopolysaccharide.



171


172

173 **Figure S2. 9-phenanthrol alleviates inflammation in experimental AP.** (A-B) Representative images of LY6G  
 174 staining and quantitative of LY6G staining; (C) Serum TNF- $\alpha$  level; (D) Serum IL-6 level. n = 6, error bars indicate  
 175 the SEM; \* P < 0.05 vs Sham or vs Control; # P < 0.05 vs Vehicle. LPS, lipopolysaccharide; TNF- $\alpha$ : tumor  
 176 necrosis factor- $\alpha$ ; IL-6: interleukin 6; AP, acute pancreatitis.




177

178 **Figure S3. *Trpm4*-knockout alleviates inflammation in experimental AP.** Representative photos of LY6G  
 179 staining and quantitative of LY6G staining. n = 6, error bars indicate the SEM; \* P < 0.05 vs Sham or vs Control;  
 180 # P < 0.05 vs Vehicle. LPS, lipopolysaccharide; AP, acute pancreatitis.




181  
 182

183 **Figure S4. *Trpm4*-plasmid overexpressed TRPM4 levels in AR42J cells.** Western blot analysis of the TRPM4  
184 expression level in AR42J cells. n = 6; error bars indicate the SEM; \* P < 0.05 vs PI-Vector. TRPM4, Transient  
185 receptor potential cation channel melastatin 4.




186  
187

188 **Figure S5. Inhibition of TRPM4 antagonizes mitochondrial dysfunction and cell death in acinar cells induced**  
 189 **by Ca<sup>2+</sup> overload. (A)** Flow cytometry analysis of Fluo-3 in mouse PACs; **(B)** Supernatant LDH levels; **(C)** Western  
 190 blot analysis of the TRPM4 expression level in mouse PACs; **(D)** Supernatant LDH levels; **(E)** ATP levels in mouse  
 191 PACs; **(F)** FRAP levels in mouse PACs; **(G)** Flow cytometry analysis of Fluo-3 in mouse PACs. n = 3-6, error bars  
 192 indicate the SEM; \* P < 0.05 vs Sham or vs Control; # P < 0.05 vs Vehicle. PACs, pancreatic acinar cells; LDH,  
 193 lactate dehydrogenase; ATP, adenosine triphosphate; FRAP: Ferric Reducing Antioxidant Power.



194  
 195

196 **Figure S6. NMDA receptor agonists exacerbate Cerulein-induced acinar cell damage.** Supernatant amylase  
197 levels. n = 6, error bars indicate the SEM. LPS, lipopolysaccharide; NMDA: N-methyl-d-aspartate.




198

199

200 **Figure S7. Inhibition of NMDAR had little effect on the expression of TRPM4 in pancreatic tissue of AP mice.**

201 **(A-B)** Western blot analysis of TRPM4 expression level in AR42J. TRPM4, Transient receptor potential cation

202 channel melastatin 4; AP, acute pancreatitis; LPS, lipopolysaccharide.



203

204

Raw western blot images

205  
206  
207  
208

FIG-1



FIG-1A-ACTIN




FIG-1A-TRPM4



FIG-1B-ACTIN




FIG-1B-TRPM4




FIG-1C-ACTIN



FIG-1C-TRPM4

209  
210  
211  
212


FIG-2



2H-ACTIN




2H-RIP3



2L-ACTIN



2L-BAX




2L-CASPASE-3

213  
214


215 FIG-3



3C-ACTIN



3C-MFN2




3C-PGC1




3C-PINK1




3I-ACTIN



3I-GRP78



3I-IRE1



3I-p-IRE1

216


217

218

FIG-4




4A-ACTIN




4A-TRPM4



4J-ACTIN




4J-RIP3




4O-ACTIN



4O-GRP78



4O-PDI



4O-P-IRE1



4O-T-IRE1

219

220

221 FIG-5




I-ACTIN




I-eIF2 $\alpha$



I-GRP78




I-IRE1



I-PDI



I-p-eIF2 $\alpha$



I-P-IRE1




K-ACTIN




K-TRPM4


222

223






230 FIG-8




8G-ACTIN




8G-PGC-1α




8H-ACTIN



8H-GRP78



8H-P-IRE1



8H-T-IRE1


231

232

233

234

FIG-S4



S4-ACTIN



S4-TRPM4


235

236

237

238

FIG-S5



S5-TRPM4




S5-β-actin

239

240

241 FIG-S7



242

243


244  
245

## Report of *trpm4* gene knockout mice

### 1. Product Information

Name	C57BL/6NCya- <i>Trpm4</i> <sup>em1</sup> /Cya
Serial Number	KOCMP-68667-Trpm4-B6N-VA
Gene	Trpm4
NCBI ID	68667
Strain	C57BL/6NCya
Type	conventional knockout

### 2. Targeting Strategy



246  
247

### 3. Genotyping strategy

Primers1: (Annealing Temperature 60.0 °C)

F1: 5'-AGTCCCTGCTCATTACTCTGGG-3'  
R1: 5'-GAGAAGGTAAAGAAGCCTGTGC-3'  
Product size: 838 bp



Primers2: (Annealing Temperature 60.0 °C)

F2: 5'-TGGGAAGTGAAGGATGAAGTCAG-3'  
R1: 5'-GAGAAGGTAAAGAAGCCTGTGC-3'  
Product size: 643 bp

Homozygotes: one band with 838 bp  
Heterozygotes: two bands with 838 bp and 643 bp  
Wildtype allele: one band with 643 bp


248  
249

#### 4. Expected Results

Genotyping	Primers1			Primers2		
	WT	Water	M	WT	Water	M
Trpm4(-/-): Homozygotes						

**Note:**

- 1) PCR was carried out in 25  $\mu$ L volume for 35 cycles under standard conditions, with all two primers listed above added to each reaction.
- 2) DNA marker: Thermo Scientific GeneRuler 100 bp DNA Ladder #SM0242



- 3) Controls used in PCR genotyping are:
  - Water control: No DNA template added.
  - Wildtype control: Mouse genomic DNA.

250

251

## 5. PCR reaction

### 5.1 DNA Extraction

#### ➤ Method One:

We recommend that using TaKaRa MiniBEST Universal Genomic DNA Extraction kit (Ver.5.0\_Code No. 9765) to gain high purity of genomic DNA.

- a. Add 180  $\mu$ L of Buffer GL, 20  $\mu$ L of Proteinase K and 10  $\mu$ L of RNase A per tail piece (2-5 mm) in a microcentrifuge tube. Be careful not to cut too much tail.
- b. Incubate the tube at 56°C overnight.
- c. Spin in microcentrifuge at 12,000 rpm for 2 minutes to remove impurities.
- d. Add 200  $\mu$ L Buffer GB and 200  $\mu$ L absolute ethyl alcohol with sufficient mixing.
- e. Place the spin Column in a collection tube. Apply the sample to the spin and centrifuge at 12,000 rpm for 2 min. Discard flow-through.
- f. Add 500  $\mu$ L Buffer WA to the spin column and centrifuge at 12,000 rpm for 1 min. Discard flow-through.
- g. Add 700  $\mu$ L Buffer WB to the spin column and centrifuge at 12,000 rpm for 1 min. Discard flow-through. (Note: Make sure the Buffer WB has been premixed with 100% ethanol. When adding Buffer WB, add to the tube wall to wash off the residual salt.)
- h. Repeat step g.
- i. Place the spin Column in a collection tube and centrifuge at 12,000 rpm for 2 min.
- j. Place the spin Column in a new 1.5ml tube. Add 50~200  $\mu$ L sterilized water or elution buffer to the center of the column membrane and let the column stand 5min. (Note: Heating sterilized water or elution buffer up to 65°C can increase the yield of elution.)
- k. To elute DNA, centrifuge the column at 12,000 rpm for 2 min. To increase the yield of DNA, add the flow-through and/or 50~200  $\mu$ L sterilized water or elution buffer to the center of the spin column membrane and let the column stand 5 min. Centrifuge at 12,000 rpm for 2 min.
- l. Quantify to genomic DNA. Eluted genomic DNA can be quantified by electrophoresis or absorbance.

#### ➤ Method Two:

A low-cost and sample method to gain rough genomic DNA.

- a. Add 100  $\mu$ L of tail digestion buffer per tail piece (2-5 mm) in a microcentrifuge tube. Be careful not to cut too much tail.
- b. Incubate the tube at 56°C overnight.
- c. Incubate the tube at 98°C for 13 minutes to denature the Proteinase K.
- d. Spin in microcentrifuge at top speed for 15 minutes. Use an aliquot of supernatant straight from the tube (1  $\mu$ L in a 12.5  $\mu$ L reaction) for PCR.

Final concentration of tail digestion buffer:

- 50 mM KCl
- 10 mM Tris-HCl (pH 9.0)

252

253

- 0.1 % Triton X-100
- 0.4 mg/mL Proteinase K

**5.2 PCR Mixture (primer concentration: 10 $\mu$ M):**

Component	x1
ddH <sub>2</sub> O	9.0 $\mu$ l
Product primer F	1.0 $\mu$ l
Product primer R	1.0 $\mu$ l
Premix Taq	12.5 $\mu$ l
DNA	1.5 $\mu$ l
Total	25 $\mu$ l

**5.3 PCR Reaction Conditions:**

Step	Temp.	Time	Cycles
Initial denaturation	94 °C	3 min	35 x
Denaturation	94 °C	30 s	
Annealing	60 °C	35 s	
Extension	72 °C	35 s	
Additional extension	72 °C	5 min	

**5.4 Relevant Reagents:**

<b>Trizma Hydrochloride Solution</b>	Sigma, Cat. No. T2663
<b>Proteinase K</b>	Merck, Cat. No. MK539480
<b>Triton X-100</b>	Sigma, T8787-50 mL
<b>Premix Taq Polymerase</b>	Vazyme, P222
<b>Agarose</b>	BIOWEST AGAROSE, REGULAR
<b>DNA Marker</b>	Thermo Scientific GeneRuler 100 bp DNA Ladder #SM0242
<b>0.5<math>\times</math>TBE</b>	Tris Bio Basic Inc, TBO194-500g
	EDTA Shanghai Sangon, 0105-500g
	Boric Acid, Shanghai Sangon, 0588-500g

254

255